



OPERATOR'S MANUAL

For the E2 Laser and Endoscopy System



OME 2000:
E2 COMPACT
MICROPROBE™ SYSTEM

TABLE OF CONTENTS

PAGE

BEGINNING

| | |
|------------------|-----|
| Warning | 1 |
| Labels | 2-4 |
| Precautions..... | 5-6 |

SYSTEM OVERVIEW

| | |
|--|----|
| System | 7 |
| Cabinet..... | 7 |
| Light Source | 7 |
| CCD Camera | 8 |
| Digital Displays and Indicators | 8 |
| Video Display..... | 8 |
| Foot Switch | 8 |
| Front Panel..... | 9 |
| Back Panel | 10 |
| Diode Laser and Principal of Operation | 11 |

GETTING READY

| | |
|---|-------|
| Site Preparation..... | 12 |
| Utilities | 12 |
| Laser Safety | 12-13 |
| Reflection Hazard | 13 |
| Tissue Protection..... | 13 |
| Explosion Hazard | 13 |
| Vapor Plume..... | 14 |
| Exposure Protection from the Aiming Beam Laser..... | 14 |
| Safe Viewing Times..... | 14 |
| Safety Features | 14-15 |

CLEANING AND STERILIZATION

| | |
|---|----|
| Cleaning the Laser Console | 16 |
| Cleaning the Laser Connector..... | 16 |
| Cleaning the Video Adapter | 16 |
| Endoscope and Probe Cleaning and Sterilization..... | 16 |

CLINICAL APPLICATIONS

| | |
|-----------------------------|----|
| Indications for Use | 17 |
| Glaucoma | 17 |
| Vitreoretinal Surgery | 17 |
| Contraindications..... | 18 |

Table of Contents (Cont.)

OPERATION

| | |
|--|----|
| Set Up and Operation..... | 19 |
| Endoscopes and Probes..... | 19 |
| Inspection of the Optical System..... | 19 |
| Eye Insertion and Videography..... | 20 |
| Eye Endophotocoagulation..... | 20 |
| Front Panel Features..... | 20 |
| Emergency Off..... | 20 |
| Standby..... | 20 |
| Enable..... | 21 |
| Aiming Beam Push Button..... | 21 |
| Laser Power Push Button..... | 21 |
| Laser Duration Push Button..... | 21 |
| Laser Output Power Control..... | 22 |
| Counter Display..... | 22 |
| Counter Reset Push Button..... | 22 |
| Preliminaries..... | 22 |
| System Turn-on..... | 23 |
| Setting of the Treatment Beam..... | 23 |
| Before Firing the Laser..... | 23 |
| Firing the Laser..... | 24 |
| Setting Single Exposure Time..... | 25 |
| Between Patient Treatment and System Turn off..... | 25 |

MAINTENANCE AND TROUBLESHOOTING

| | |
|--|-------|
| Endoscope Maintenance..... | 26 |
| Laser Maintenance..... | 26 |
| Laser Power Calibration..... | 27 |
| Maintenance and Troubleshooting Guide..... | 28-30 |
| Operator Replaceable Parts..... | 30 |

TECHNICAL SPECIFICATIONS..... 31-32

EMC GUIDELINES.....33-36

BIBLIOGRAPHY.....37-42

ECP TREATMENT SUGGESTIONS..... .43

QUICK SET UP GUIDE..... .44

E2 COMPACT MICROPROBE™

MICROENDOSCOPE SYSTEM

Endo Optiks, Inc.
39 Sycamore Avenue
Little Silver, NJ 07739-1208
USA

Tel: 001 732-530-6762 Fax: 001 732-530-5344
E-mail: info@endooptiks.com
Website: <http://www.endooptiks.com>



Advena Ltd. Pure Offices, Plato Close, Warwick, CV34 6WE UK.

WARNING: NO UNAUTHORIZED USE OF LASER. The user of the E2 MicroProbe™ should be thoroughly trained in the applicable procedure. Furthermore, failure to read and thoroughly understand the content of this Operators Manual may result in serious injury to the patient or user. It is essential to follow the instructions contained in this manual which pertain to the E2 MicroProbe™ and accessories used in conjunction with the procedures. Failure to follow these instructions may result in damage to the E2 MicroProbe™ or malfunction of the E2 MicroProbe™.

CAUTION: Endo Optiks restricts the sale of the E2 MicroProbe™ to a physician or on order of a physician.

BEGINNING

Labels

The following labels are affixed to the E2 Microprobe™ system. The title, the part number, and the location on the E2 Microprobe™ are given for each label.

Label

Location



Identification

P/N 3840406
Rear Panel, near the top right corner.



Meaning: Type BF equipment

Protection Against Electric Shock

P/N L1012
Incorporated into the Identification label

Supply Rating

115/240 Volts, 60/50Hz, 6.0/3.0A
Use 2 Type T6.3A, 250V

Fuse Replacement Label

P/N 0711111
Rear Panel, directly below AC receptacle



Meaning: Dangerous Voltage

P/N L1007
Incorporated into the Identification label

BEGINNING



LASER APERTURE

P/N L1003C

Front Panel



Caution: Consult Accompanying Documents

P/N L1009

Incorporated into the Identification label



LASER WARNING

P/N L1002C

Top, Front

This product conforms to the applicable Requirements of US 21 CFR, subchapter J, FDA Laser Notice 45, Laser Guide: 1995, and IEC 60825-1:2007

CONFORMITY & APPLICABLE STANDARDS

P/N L1016


Rear Panel



LASER STOP

P/N L1017


Front Panel



Conforms to European Medical Device Directive 93/42/EC.

P/N L1003

Incorporated into the Identification label



Safety Agency Approvals (MET MARK)

P/N E112374

Rear panel, bottom left, under AC receptacle



Caution: HOT

P/N: L1009A

Front Panel, next to light connector when 300W light installed



This symbol has been attached to the equipment or, in the case that this is not possible, on the packaging, instruction literature and/or the guarantee sheet. By using this symbol it states that the device has been marketed after August 13th 2005, and implies that you must separate all of its components when possible, and dispose of them in accordance with local waste disposal legislations.

- Because of the substances present in the equipment, an improper use or disposal of the refuse can cause damage to human health and to the environment.
- With reference to RAEE (Registry of Electrical and Electronic Apparatuses), it is compulsory not dispose of the equipment with normal urban refuse, arrangements should be instigated for separate collection and disposal.
- For more detailed information about recycling of RAEE, please contact your local waste collection body.
- In case of illicit disposal, sanctions will be levied on transgressors.

To achieve proper grounding reliability, a power supply plug must be fully inserted into a receptacle marked "HOSPITAL GRADE"

Hospital Grade Power Cord

P/N: Hospital Grade Labels

Above Power Receptacle On Rear Panel

BEGINNING

Precautions

To prevent fire or shock hazard, do not expose the unit to rain or moisture.

Dangerously high voltages are present inside the E2 Microprobe™. Do not open the cabinet. Refer servicing to qualified personnel only.

In the event of a malfunction or when maintenance is necessary, consult:

Endo Optiks

39 Sycamore Ave., Little Silver, NJ, USA.

Tel: 001 732 530 6762, Fax: 001 732 530 5344

Email: info@endooptiks.com

On safety

- Operate the unit on the designated AC voltage only.
- The Fuse Replacement Label indicates operating voltage and is located adjacent to the mains fuse holder in the rear of the cabinet.
- Should any solid object or liquid fall in, unplug the unit and have it checked by qualified personnel before operating it any further.
- To disconnect the AC power cord, pull it out by grasping the plug, never pull the cord itself.
- The outlet shall be installed near the equipment and shall be easily accessible.

Warning

This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the E2 Microprobe or shielding the location.

Endoscopes and Probes

This device is intended to be used in conjunction with Endo Optiks endoscopes and probes ONLY and to assure safety should not be connected or used with any other devices.

Disconnect the endoscope from the system by grasping the connectors only. Do not grasp or bend the endoscope jacketing for this may break the glass fibers that are enclosed within the jacketing.

CAUTION: The heat generated by the light source varies and can be hot causing the light connector, on the endoscope, to become hot. To prevent injury to the operator the light intensity should be turned down so the connector can cool before removing it from the system.

On installation

- The E2 MicroProbe™ should be used in a Hospital or Clinical setting only. It should be used indoors only under the environmental conditions stated on Page 36 of this document.
- Allow adequate air circulation to prevent internal heat build-up.
- Do not place the unit on surfaces (rugs, blankets, etc.) or near materials (curtains, draperies) that may block the ventilation holes.
- Do not install the unit in a location near heat sources such as radiators or air ducts, or in a place subject to direct sunlight, excessive dust, mechanical vibration or shock.

On cleaning

To keep the unit looking brand-new, periodically clean it with a mild detergent solution. Never use strong solvents such as thinner or benzine, or abrasive cleansers since they will damage the cabinet. As a safety precaution, unplug the unit before cleaning it.

On sterilization before use

WARNING - The endoscopes and probes must be sterilized before use. Please refer to the instructions provided with each device.

SYSTEM OVERVIEW

System

The Endo Optiks E2 MicroProbe™ is the principal component in a new portable laser and endoscopy system. The complete system consists of the therapeutic laser, the endoscope, the monitor and the footswitch. This compact unit creates the opportunity to simultaneously image and photocoagulate the ciliary processes through a corneal incision. It is especially indicated for the safe and effective treatment of glaucoma in combination with cataract surgery. Important vitreo-retinal applications can be realized. It can be used for the contact and non-contact excision, hemostatis, incision and vaporization of soft tissue.

Cabinet

The compact laser and endoscopy cabinet houses a xenon light source, a therapeutic laser and a CCD camera. The laser output, pulse width, light and aiming beam intensity are controllable from the Front Panel. The parameters are displayed on Front Panel digital displays and lighted status indicators. For safety there is an emergency shutoff button. The Rear Panel features connectors to any video monitor, VCR or video printer. A Foot Pedal enables hands-free operation.

Light Source

The xenon light source is used to provide light to the endoscope. The intensity of light can be adjusted from the Front Panel, or with the foot switch.

CAUTION: The heat generated by the light source varies and can be hot causing the light connector, on the endoscope, to become hot. To prevent injury to the operator the light intensity should be turned down so the connector can cool before removing it from the system.

SYSTEM OVERVIEW

CCD Camera

The CCD camera is used to process the image obtained by the fiberoptic endoscope and display it on the video display. There is a Video Camera BNC connector Input and Outputs located at the Back Panel. The Video Camera Cable Output can be plugged into the Video Camera Cable Input or a remote Video Camera Cable Input can be used.

Digital Displays and Indicators

The Display consists of four Light Emitting Diodes and four Lighted Indicators. They are as follows:

| <u>FUNCTION</u> | <u>DISPLAY TYPE</u> | <u>Key*</u> |
|------------------------|----------------------------|--------------------|
| Laser Power | Digital | 1 |
| Laser Active | Lighted | 2 |
| Laser Duration | Digital | 3 |
| Laser Standby | Lighted | 4 |
| Laser Enable | Lighted | 5 |
| Laser Aiming Power | Digital | 6 |
| Laser Cool Down | Lighted | 7 |
| Laser Shot Counter | Digital | 8 |

* Please refer to the Diagram labeled Front Panel on Page 8.

Video Display

The video display can be any high resolution monitor such as the Sony LMD-1530MD and is used for displaying the endoscopic image. The video outputs are located at the Back Panel and can be utilized for recording the endoscopic image onto any video recording format such as the U.S. standard NTSC or the European standard PAL. There is an S-Video (Y/C Out) and 4 Video Out Connectors). All are BNC connectors (75 ohms terminated).

Foot Switch


The footswitch is used to activate the laser. Some footswitch models can also be used to vary the illumination intensity of the xenon light source. The footswitch connector is located at the Back Panel.

SYSTEM OVERVIEW

Front Panel

The Front Panel contains the color coded switches, digital displays and indicators used to control the functioning and show the status of the E2 Microprobe™. The functions and status indicators are:

| <u>KEY</u> | <u>COLOR</u> | <u>FUNCTION</u> |
|-------------------|---------------------|--|
| 1 | | Laser Power Digital Display |
| 2 | | Laser Active Lighted Indicator |
| 3 | | Laser Duration Digital Display |
| 4 | | Laser Standby Lighted Indicator |
| 5 | | Laser Enable Lighted Indicator |
| 6 | | Laser Aiming Power Digital Display Counter |
| 7 | | Laser Cool Down Lighted Indicator |
| 8 | | Laser Shot Counter Digital Display |
| 9 | | Laser Output Connector |
| 10 | Yellow | Laser Power Up and Down Switch |
| 11 | Blue | Laser Duration Up and Down Switch |
| 12 | Green | Laser Enable Switch |
| 13 | Red | Aiming Beam Intensity Up and Down Switch |
| 14 | White | Illumination Intensity Up and Down Switch |
| 15 | | Xenon Lamp Out |
| 16 | | Counter Reset |
| 17 | Red | STOP (Big Red Button) |
| 18 | | Camera |



Front Panel Features

SYSTEM OVERVIEW


Back Panel

The Back Panel contains the On -- Off Key Switch, the Mains Power Inlet, the Fuses, the Foot Switch Connector, the Video Connectors, the Remote Interlock and the Remote Communications Port. The locations are:

KEY

FUNCTION

| | |
|----|---|
| 1 | Video Camera Cable (located under the cabinet foot) |
| 2 | Video Camera Cable Input |
| 3 | S-Video (Y/C Out) |
| 4 | Foot Switch |
| 5 | Remote Interlock |
| 6 | Laser On -- Off Key Switch |
| 7 | RS-232 Remote Communications Port |
| 8 | ON/OFF Switch and Mains Power Inlet and Fuses |
| 9 | Video Out |
| 10 | Video Out |
| 11 | Fuse for Lamp Power Supply |
| 12 | Video Out |
| 13 | Video Out |



Back Panel Layout

SYSTEM OVERVIEW

Diode Laser and Principal of Operation

Diode lasers are small semiconductor devices consisting of a sandwich of gallium-aluminum arsenate crystalline materials and end mirrors. The electrons of the crystal are raised to an excited energy state by an electrical current. When the electrons return to their original energy state, photons are emitted. In the laser cavity, these emitted photons are trapped between the cladding layers and the end mirrors. When a photon passes close to an excited electron, the electron will be stimulated to emit another photon that is identical in wavelength and phase to the first. This amplification process continues, increasing the number of active photons as the photon light beam is reflected back and forth between the cavity mirrors. One of these mirrors releases a percentage of the energy hitting its surface resulting in the infrared laser light.

The frequency of the diode laser employed in this system is 810nm. This laser is intended only for the use of physicians who are trained in operation of laser photocoagulators. Training in the therapeutic use of lasers is available in medical courses and seminars offered throughout the world.

GETTING READY

Site Preparation

The E2 Microprobe™ system has no special electrical or water requirements. Standard wall voltage can be utilized and no special cooling is required because a solid state diode laser is used.

Utilities

Electrical: The E2 Microprobe™ can be configured to operate from a power source in the range of either 115 volts AC or 240 volts AC, 60/50Hz, 690/720 watts maximum. A standard grounded AC outlet is sufficient. Verify that the voltage indicated on the fuse holder, below the A/C receptacle, at the back of the laser console matches the actual line voltage before the instrument is plugged in.

WARNING

This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the E2 Microprobe or shielding the location. Refer to EMC Guidance in the Technical Specifications section.

WARNING

No modification of this equipment is allowed.

WARNING

Do not modify this equipment without authorization of the manufacturer.

WARNING

If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.

Laser Safety

WARNING - Never look directly into the laser aperture (output port) or fiberoptic when power is applied. Severe eye damage could occur.

As a precaution against accidental exposure to the output beam or its reflection, all persons in the vicinity during operation of the photocoagulator must wear laser safety glasses. The only exception is the surgeon if he or she is looking through a delivery system which is protected by an internal laser filter.

Never look directly into the laser beam or its reflection, even with the laser safety glasses in place. The laser safety glasses must give protection at 810NM.

Note: The Endo Optiks diode laser safety glasses have been tested and are compliant with EN60207-EN60208 standards. The 810nm wavelength and Optical Density (O.D.) of 5+ are clearly indicated on the laser safety glasses. Safety glasses that are obtained from any vendor other than Endo Optiks should be compliant with these standards. Some multiwavelength glasses will block the aiming beam and laser firing indicator, on the front panel of the system and should not be used.

Reflection Hazard

Smooth objects will reflect a diode laser beam. Reflection hazards can exist beyond the laser beam aperture. Regardless of the surface, reflection is a potential hazard when the laser strikes a nonabsorbent surface such as a metallic instrument. There are two types of reflections which can occur: Specular and diffuse.

Specular reflections: Specular reflection occurs when the size of the surface irregularities is less than the wavelength of the incident radiation. In this instance the laser beam reflects at the same angle as the angle of incidence. The danger occurs when a concentrated beam of laser light could accidentally strike unintended tissue or the eyes causing significant thermal damage.

Diffuse reflections: Diffuse reflection occurs when surface irregularities are randomly oriented and are much greater than the wavelength of the incident light. In this instance the reflected laser beam is disorganized and scatters in many directions though a laser light may still strike an unintended target surface. But the power density will be reduced to a degree that will cause minimal, if any, thermal effect.

Tissue Protection

WARNING - Never place your hand or other objects in the pathway of the laser beam. Severe burns could occur.

Explosion Hazard

WARNING - Never operate the laser in the presence of flammable anesthetics. An explosion and fire could occur.

Do not use the laser in the presence of flammables or explosives such as volatile anesthetics, alcohol, certain surgical preparation solutions and other such substances.

The E2 MicroProbe™ laser should be disabled whenever the system is not actually in use during treatment of a patient. To safely disable the laser the key should be removed from the rear panel of the system. This is to prevent accidental laser discharge if the Foot Pedal is inadvertently depressed

GETTING READY

Vapor Plume

CAUTION: Laser plume may contain viable tissue particulates; however the endoscope is fired in a closed environment in both the anterior and posterior segments of the eye so the plume does not escape into the atmosphere. There are no possible hazards associated with the laser plume.

Exposure Protection from the Aiming Beam Laser

CAUTION: The maximum safe power that can be visualized indefinitely by the retina is 0.39 microwatts (Class 1). The diode laser that produces the red aiming beam is a Class 3R laser, producing up to 1,500 microwatts. Direct visualization of Class 3R laser radiation by a patient during treatment has not been shown to be safe. To protect the patient from possible retinal damage use the lowest practical aiming beam intensity during treatment. The diode laser produces a red diode aiming beam with power variable from barely visible to 1.5 milliwatts maximum. The safe (Class 1) exposure duration limit at maximum power level of 1.5 milliwatts is 0.047 seconds.

Safe Viewing Times

| Class | Description | Type |
|-------|--|------------------------|
| 1 | Safe under conditions of operation. | |
| 2 | Eye protection by aversion response (blink). | |
| 3 | Minimum safe viewing time. | Aiming Laser |
| 4 | Use requires extreme caution. | Treatment Laser |

Safety Features

Key Lock Switch: The laser can be turned on only with the proper key to close the master switch. The key cannot be removed in the On position and the laser will operate only when the key is in place. When the Key Switch is rotated to the On position, power is available to the instrument. For safety purposes, the key should be removed when the laser is not being used.

Enable Function: Even when the E2 MicroProbe™ is activated by turning the key on, the Foot Pedal will not fire the laser. The laser cannot be fired unless the ENABLE control is pressed on. This control takes the laser out of the

STANDBY mode and allows laser firing by footswitch activation. This safety feature requires a deliberate second action by the Operator of the laser prior to its being fired.

Protective Housing: The diode laser has a protective housing which prevents unattended human access to laser radiation above safe limits. The user shall not perform maintenance or service operations inside the protective covers. This housing is to be opened only by an Endo Optiks Service Representatives or an Endo Optiks certified technician.

Protective Jacketing: The laser fiber inside the endoscope is protected by jacketing. In the event that the laser fiber breaks, the jacket protects the user from laser radiation emission.

Viewing Optics: When employing the laser endoscope in the Video mode, direct viewing of the laser and ocular tissues is not an issue, because the laser energy cannot be emitted through the video display.

Room Door Interlock: The E2 MicroProbe™ will be switched into standby mode when the treatment room door is opened during use.

System Labels: Appropriate warning labels have been mounted in specified locations on the instrument to indicate conditions under which the user could be subjected to laser radiation.

STOP: Depressing the STOP (the Big Red Button) will immediately disable the laser and reset all the treatment parameters to zero. The E2 MicroProbe™ will be in the stand-by mode with the power on.

WARNING

It is not intended that the user perform any maintenance or service operations inside the protective housing of the diode laser. All internal adjustments, cleaning, maintenance, and servicing shall be performed only by an Endo Optiks Service Representative or an Endo Optiks certified technician.

CAUTION

Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

CLEANING AND STERILIZATION

Cleaning the Laser Console

To clean the external surface of the Laser Console, wipe using a cloth dampened with a noncaustic cleaning solution such as soap and water, isopropyl alcohol or a "hospital grade" disinfectant. Do not spray or put cleaning agents directly on the system. Dry with a clean, dry cloth or allow to air dry.

Cleaning the Laser Connector

Before connecting an endoscope to the system, clean the tip of the Laser Connector with alcohol and lens cleaning tissue. Moisten the lens tissue with alcohol and gently wipe it once across the tip of the Laser Connector.

Cleaning the Video Adapter

Before connecting an endoscope to the system, clean the tip of the Video Connector with alcohol and lens cleaning tissue. Moisten the lens tissue with alcohol and gently wipe it once across the tip of the Laser Connector. This will prolong the life of the endoscope and maximize image quality.

Clean the Video Adapter with a lens cleaning tissue. Do not use alcohol for this may harm the laser filter that is built into the Video Adapter.

Endoscope and Probe Cleaning and Sterilization

These devices are delivered non-sterile and are designed to be reusable in accordance with the constraints described in the information supplied. The Endoscope should be sterilized before each use. If the endoscope is sterilized but not cleaned after each use, you will almost certainly destroy the tip of the endoscope the next time you fire the laser.

FULL INSTRUCTIONS FOR CLEANING AND STERILIZATION ARE PROVIDED WITH EACH ENDOSCOPE OR PROBE.

CLINICAL APPLICATIONS

WARNING

The MicroProbe™ system is intended solely for use by physicians trained in the operation of laser photocoagulators.

Indication For Use

The ophthalmic laser endoscope is indicated for the Endo intraoperative photocoagulation of the ciliary processes in the treatment of glaucoma, proliferative retinopathies, retinal detachment, focal treatment of retina or choroidal disorders, and for evaluation of the internal ocular structures in patients with dense opacifications of the anterior segment which do allow a posterior view.

Glaucoma

This instrument is indicated for the treatment of glaucoma in patients who have failed with conventional topical and systemic medications, or previous laser photocoagulation, or trabeculectomy and other filtering procedures, or cyclocryotherapy or other cyclodestructive procedures.

The endophotocoagulation of ciliary processes under direct endoscopic view is highly controllable, i.e. titratable, and has been demonstrated to be effective in the treatment of glaucoma.

Vitreoretinal Surgery

The endoscopically controlled endophotocoagulation that is possible with the ophthalmic laser endoscope is also useful for endophotocoagulation:

- During vitreous surgery to produce chorioretinal scar around retinal breaks or retinotomy sites
 - To perform intraoperative panretinal photocoagulation (PRP) in proliferative retinopathies
 - To perform intraoperative retinal photocoagulation on a scleral buckle
 - To perform intraoperative photocoagulation around focal neovascularization

CLINICAL APPLICATIONS

Contraindications

CAUTION: It is NOT advisable to use the E2 Microprobe if the following contraindications are present.

USE CONTRAINDICATED IN:

- An albino eye

WARNING

The diode laser and endoscopic photocoagulation are intended solely for use by physicians trained in the operation of photocoagulators.

WARNING

Excessive number or short pulse durations may result in choroidal hemorrhage.

WARNING

Excessive laser power settings may result in retinal holes and retinal hemorrhages.

WARNING

Photocoagulation of ciliary processes may result in anterior or posterior segment hemorrhage, uveitis or hypotony. A diode laser has an effect on tissue similar to that of the krypton laser. The longer wavelengths of krypton and diode lasers relative to argon lasers may result in the following situations:

- An increased uptake by the choroid which may cause
 - a. Increased patient discomfort
 - b. Increased choroidal edema
 - c. Increased choroidal hemorrhage
 - d. Reduced accommodation
- Poor uptake by hemoglobin in management and treatment of retinal capillary aneurysms, retinal arteriolar macroaneurysms and retinal capillary hemangioma.
- Poor uptake in a hypopigmented eye.

OPERATION

Set Up and Operation

The E2 MicroProbe™ system must be set up and operated with care. The endoscopic laser probe is a fragile component and great care should be used when handling it.

Training on the safe operation of the system should be provided at the time of installation. This includes complete in-service training. This training session is usually conducted in one day. From time to time it may become necessary to provide retraining. This can be provided by contacting your sales representative.

Endoscopes and Probes

Before each use, the endoscope or probe should be thoroughly inspected in order to detect any irregularity or abnormality. If these are noted, the endoscope or probe should not to be used and a new one should be employed.

Inspection of the Optical System

The endoscope allows two functions. The first is optical viewing and the second is delivery of therapeutic laser energy. Both of these functions can be determined prior to insertion into the eye. It is essential that the viability of the fiberoptic pathways be demonstrated prior to use and this may be done in a rather simple manner.

To check the integrity of the optical viewing system, direct the distal end of the endoscope to a point of reference, such as the focus ring at the video adapter for the image fiberguide on the Front Panel of the E2 MicroProbe™. A clear view should be obtained.

The aiming beam provides a good method of checking delivery of the therapeutic laser energy since it passes down the same delivery system as the therapeutic laser beam. If the aiming beam spot is not present at the distal end of the delivery system, its intensity is reduced or it looks diffused, this is a possible indication of a damaged or not properly working delivery system.

OPERATION

Eye Insertion and Videography

Insertion: The E2 MicroProbe™ endoscope may be inserted in the eye through the pars plana through a standard vitrectomy incision or may be inserted through a limbal incision that has had previous or concurrent surgery.

Observation: Observation of the intraocular structures occurs by viewing of the high resolution video monitor. The internal structures of the eye from the posterior aspect of the iris, ciliary body, pars plana, peripheral retina and more posterior retina may be imaged.

Videography: Video outputs are included in the E2 MicroProbe™ system. Simply attach a video recorder to these outputs, insert a video tape and depress the PLAY & RECORD buttons simultaneously. The video image will be recorded. **Digital video recording devices can also be utilized.**

Eye Endophotocoagulation

Endophotocoagulation of the ciliary processes, pars plana, peripheral retina and more posterior retina is possible through this instrument.

Front Panel Features

The following controls and indicators are located on the Front Panel:

Emergency Off (STOP Button)

Depressing the Emergency Off (Big Red Button) de-energizes the lasers. All laser parameters are then reset to zero. The power supply remains on. The switch does **NOT** de-energize the entire unit. It is not intended as a normal shutdown control.

Standby

The E2 MicroProbe™ is in the Standby mode when first turned on. The Foot Pedal is disabled and therefore no laser treatment beam is possible. The aiming beam, however, is available. The Standby lighted status indicator will illuminate to indicate that the system is in the Standby mode.

OPERATION

WARNING

Except during the actual treatment, the system should always be in the Disabled mode. This can be achieved by removing the key from the rear panel. Maintaining the system in the Disabled mode prevents accidental exposure if the Foot Pedal is inadvertently depressed.

Enable

Key # 12, page 8

The E2 MicroProbe™ is put in the Enable mode by pressing the Enable switch. In Enable mode, the Foot Pedal is activated; depressing the footswitch will deliver the treatment beam to the target. The Enable indicator will illuminate to alert the operator that the system is in the Enable mode.

Aiming Beam Push Button

Key # 13, page 8

The Aiming Beam up and down push buttons are used for controlling the aiming beam intensity. Press the up button for higher intensity and the down button for lower intensity.

Laser Power Push Button

Key # 10, page 8

The Laser Power up and down push buttons are used for selecting laser output power. Press the up button for higher power and the down button for lower power. The power changes in 10 milliwatt increments.

Laser Duration Push Buttons

Key # 11, page 8

The Laser Duration push buttons are used for selecting laser exposure time. Exposure time is the amount of time the laser produced a treatment beam and is measured in seconds. Preset exposure times range from 0.05 to 2.0 seconds; and, continuous wave is noted as “Con”. Press the Laser Duration up button to increase the exposure time; press the Laser Duration down button to decrease the exposure time. When the time indicated on the setting has elapsed, the treatment beam stops even if the Footpedal is still depressed (unless a continuous laser exposure has been selected).

OPERATION

Laser Output Power Control

The initial laser output power is permanently recorded by the E2 Microprobe™. If the subsequent laser output power varies by 20% or more, the Laser Output Power Meter will start blinking.

Counter Display

Key # 8, page 8

The Counter display indicates the number of exposures delivered to the target tissue since the last reset.

Counter Reset Push Button

Key # 16, page 8

Press the Reset button on the front control panel to zero the counter display.

E2 MicroProbe™ Preliminaries

- Ensure the E2 MicroProbe™ is connected to the mains.
- Ensure the Footswitch is properly connected.
(See Connection Points diagram on page 1-1, Topic 1.1 System)
- Ensure the Laser Endoscope is connected to the Front Panel.
 1. Attach the Laser Connector to the Laser Aperture receptacle on the Front Panel. Screw the connector on.
 2. Attach the Video Connector to the Video Adapter on the Front Panel.
 3. Attach the ACMI Light Connector to the Light Source on the Front Panel.

The E2 MicroProbe™ is now ready for Turn-on.

WARNING

To prevent explosion and fire do not operate the laser in the presence of flammables or explosives such as volatile anesthetics, alcohol or other such substances. Do not look into the Laser Output at any time. Severe eye damage could occur.

E2 MicroProbe™ System Turn-on

1. Turn the mains power switch “On”. A system Self-Test begins.
2. The audible Pre-Laser Firing (PLF) Warning Alarm will remain ON for approximately 10 seconds.
3. After the PLF Warning Alarm ceases, the system will go into the Standby mode.
4. The STANDBY light on the display panel illuminates after the PLF Warning Alarm and the COOLDOWN light go off (if it illuminates).

The E2 MicroProbe™ is now ready for Setting of the Treatment Beam.

E2 MicroProbe™ Setting of the Treatment Beam

To set the treatment beam:

1. Insert the key in the Key Switch and turn to the "On" position.
2. Depress the Laser Power push button until the desired laser power appears in the Laser Power display. The lowest practical setting should always be used.
3. Depress the Laser Duration push button until the desired exposure time appears in the Laser Duration display.
4. Press the Enable push button and the green LED status indicator will light. The 810 NM laser will not fire until the Foot Pedal is used.

SAFETY PRECAUTION

The LASER DURATION and the LASER POWER must be set to get an aiming beam.

Before Firing the Laser

1. Adjust the Light Source push buttons to obtain a satisfactory video image. If the Light Source is set too high, the aiming beam will appear diminished.
2. Bring the target video image into clear focus by rotating the star wheel on the Video Adapter while viewing the target through the endoscope.
3. Place the tip of the endoscope near a white paper. Press the Aiming Beam Intensity push button up and verify that the red aiming beam intensity increases on the paper and the monitor. An audible modulating tone indicates an aiming beam firing.

OPERATION

Firing the Laser

1. Target the tissue to be treated with the aiming beam. The lowest practical setting should always be used.

CAUTION

The maximum safe power that can be visualized indefinitely by the retina is 0.39 microwatts (Class 1). The diode laser that produces the aiming beam is a Class 3R laser. Direct visualization of Class 3R laser radiation by a patient during treatment has not been shown to be safe. To protect the patient from possible retinal damage, use the lowest practical aiming beam intensity during treatment. Start at zero and go up to the minimum amount that is practical. Endo Optiks recommends an aiming beam setting of 40.

2. Depress the Foot Switch to deliver the treatment beam to the tissue.

SAFETY PRECAUTION

The endoscope must be connected to get an aiming beam and to fire the laser.

NOTE

Early release of the Foot Switch will terminate the treatment beam before the indicated exposure time has elapsed. This will alter the therapeutic effect. To ensure that the treatment beam has exposed for the proper duration, hold the Foot Switch down until the selected exposure time has expired. The treatment beam exposure can be halted by releasing the Foot Switch before the end of the selected exposure time.

3. If surgery is interrupted, remove the key from the rear panel to disable the laser.

WARNING

Do not leave the laser unattended when it is on. Except during actual treatment, the system must always have the key removed from the rear panel. Maintaining the system with the key removed prevents accidental laser exposure if the Foot Switch is inadvertently depressed.

OPERATION

Setting Single Exposure Time

To deliver a single exposure to the target tissue:

1. Using the Laser Duration push button, adjust the exposure time to the desired setting. The selected exposure time will appear in the Laser Duration display.
2. Select the Enable mode with the Enable push button as indicated by the Enable status indicator. Each depression of the Foot Switch will deliver only one exposure at the selected exposure time. The yellow display will count down the exposure time.

Between Patient Treatments and System Turn Off

When treatment is complete or is to be interrupted:

1. Remove the key from the rear panel.
2. If several patients are to be treated, always return all Front Panel controls to their minimum settings. Disinfection and sterilization procedures are necessary and required between each surgical operation. The E2 Microprobe™ should be de-energized for cleaning.
3. When all treatments are concluded, turn the Key switch to the “O” (off) position. Remove the key to prevent unauthorized use of the laser system.

MAINTENANCE AND TROUBLESHOOTING

Endoscope Maintenance

Handling: The endoscope fiberoptic cable can be damaged or fractured if mistreated. Be careful not to step on, kink, tightly coil, pull the fiber, or catch it on any equipment. Do not clamp the fiber with a hemostat. Such stress can damage the fiber which may result in accidental laser exposure to the operating personnel or the patient. Undue physical stresses should not be applied to this instrument as it will shorten its useful life. The endoscope should be inspected monthly for physical damage to ensure that the endoscope is within the compliant criteria stated within the manual.

Sterilization: Sterilization methods for the endoscopes and probes are described in the inserts provided with each device.

Replacement: The endoscope may be reused until the optical component or the laser delivery component proves to be inadequate, as previously described. If either of these two variables proves to be inadequate, the endoscope should be discarded (dispose of in accordance with medical waste procedures) and a new one employed.

Laser Maintenance

LASER EMERGENCY TURN OFF

To de-energize the laser system in an emergency situation: Depress the STOP push button.

Introduction: The E2 MicroProbe™ system has been designed to give lifelong, trouble-free service. Service can be obtained by shipping the unit to our facility for system performance and laser safety checks.

User Maintenance for the Laser: There are no maintenance procedures to be performed by the user to ensure the intended performance of the laser, other than measuring the power of the laser.

MAINTENANCE AND TROUBLESHOOTING

Laser Power Calibration: This system is calibrated by the manufacturer and cannot be calibrated by the user. Opening the enclosure will void the warranty on this product. The user can measure the power of the laser to ensure the system is operating correctly. The laser power should be checked on an annual basis. Laser power is measured using a laser power meter capable of measuring an 810nm wavelength. Using a meter at any other wavelength will produce inaccurate results. The power should be measured at the tip of a clean endoscope. The tip should be about ¼” from the meter receiver. The power reading should be +/-20% of what the system is set for. The power output reading at the tip of the microendoscope is calibrated to be ± 20% of the power setting. If the laser power reading is more or less than 20% of the power setting the laser power display will flash and Endo Optiks should be called to schedule servicing.

NOTE: Insufficient laser power could be a result of debris on the tip of the microendoscope. Particular care should be employed at all times to ensure that the tip is clean and free of tissue or any other debris. Reduced laser power could also be caused by degradation of the laser fiber in the microendoscope and a new microendoscope should be employed.

WARNING

LASER SERVICE is to be done ONLY by Endo Optiks service engineers. Possession of service instructions or tools does not authorize repair or modification of this instrument.

WARNING

As a precaution against accidental exposure to a laser output beam or its reflection, anyone servicing the E2 Microprobe™ should wear specially tinted laser safety glasses.

WARNING

If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.

MAINTENANCE AND TROUBLESHOOTING

E2 MicroProbe Maintenance

The ONLY user replaceable maintenance items are the mains fuses. They are located on the Back Panel of the E2 Microprobe™. The mains fuse holders are in the lower left hand corner.

For 115 volts

Supply Rating

115 Volts, 50/60~, 6A

Use 2 type T6.3A, 250V

Schurter P/N FSM 034.2521

For 240 volts

Supply Rating

240 Volts, 50/60~, 3.0A

Use 2 type T6.3A, 250V

Schurter P/N FSM 034.2521

Troubleshooting

The treatment laser has been subjected to rigorous mechanical and environmental tests. In the unlikely event that the laser fails to operate properly, the Troubleshooting Guide will help you locate and correct the malfunction. Should a major malfunction occur a service representative must be contacted.

Before proceeding to the Troubleshooting Guide, check the following items:

1. Electrical Power.

- Verify that the main electrical power on the wall is on.

2. Main Power Circuit Breaker.

- Verify that the circuit breaker is in the On Position.

MAINTENANCE AND TROUBLESHOOTING

E2 Troubleshooting Guide

| DIFFICULTY | PROBABLE CAUSE | SUGGESTION |
|---|--|---|
| System Does Not Turn On | Laser Fuse Is Blown | Replace Fuse |
| | Main Power Circuit Breaker Is In The Off Position | Reset The Main Power Circuit Breaker To The On Position |
| Inadequate Aiming Beam | Aiming Beam On Low Setting | Increase Aiming Beam Setting |
| No Aiming Beam | Laser Duration and Power Not Set | Laser Duration and Power Must Be Set Before the Aiming Beam Can Be Enabled |
| | Endoscope Not Connected Properly | Check Connection of Endoscope |
| | Emergency Stop is Depressed | Release Emergency Stop Button |
| No Laser Treatment Beam | Laser Enable Mode Not Selected | Depress Enable Push Button |
| | Emergency Stop is Depressed | Release Emergency Stop Button |
| No Laser Treatment Beam or White Light or Video | Delivery System Not Properly Plugged In | Check Endoscope Connections |
| Inadequate Laser Treatment Power | Low Level Of Laser Transmission | Use New Endoscope |
| | Internal Problem | Call Endo Optiks Service |
| Blurry Video Image | Video Adapter Out of Focus Tissue on Endoscope Tip Video Connector Not Seated In Adapter Or O-Ring Missing | Adjust Video Adapter Focus Clean Endoscope With Merocel Wipe Check Connector Seating And O-Ring |
| Black spot(s) on Image | Dirty Video Adapter Lens | Clean Video Adapter Lens |
| | Dirty Endoscope Video Connector | Carefully Clean Center of Endoscope Video Connector |
| No Video Image | Monitor is turned Off | Turn Monitor On |
| | Video Cables Not Properly Connected. | Check For Proper Cable Attachment. |
| | Signal Not Selected | Depress A or B on Monitor |
| Cool Down Light Stays On | Check For Proper Airflow Under Console | Turn Off Unit, Clear Obstruction and Allow Unit to Cool. Restart Unit. |
| | Emergency Stop Switch Pressed In Key Is Not In "1" Position | Turn Emergency Stop Switch Clockwise ¼ Turn Turn Key To "1" Position |
| | Endoscope Is Not Connected | Connect Endoscope To System |
| | Software Malfunction | Turn Off Unit. Disconnect and Reconnect Endoscope. Restart Unit. |

MAINTENANCE AND TROUBLESHOOTING

E2 Troubleshooting Guide (cont.)

| DIFFICULTY | PROBABLE CAUSE | SUGGESTION |
|---|--|--|
| 777 Showing on Laser Power Display | Not Enough Amperage Available From Power Source Front Panel Button Sticking Footswitch Has An Electrical Short | Plug Unit Into a Non-overloaded Outlet Press Each Button To Verify It Clicks Twice Footswitch Needs to Be Rewired – Call Endo Optiks Service |
| No White Light or Intensity Not Adjustable | Footswitch Inadvertently Depressed. White Light Source Malfunction | Reposition Foot Call Endo Optiks Service |
| Laser Power Indicator Flashes “FFF” Or Goes Out When Laser Is Fired | Footswitch Cable Is Twisted Where It Enters The Footswitch Internal Switch In Footswitch is Defective The footswitch is detected as depressed on startup | Recycle the AC power Call Endo Optiks Service Call Endo Optiks Service Remove foot from footswitch when powering up system. |
| All displays show all “P” | The measured laser power is not within the expected limits | Recycle the AC power Call Endo Optiks Service |
| All displays show all “C” | The current supplied to the laser is not within the expected limits | Recycle the AC power Call Endo Optiks Service |
| All displays show all “S” or “5” | The Main Microcontroller does not receive the proper acknowledgement from the Slave Microcontroller | Recycle the AC power Call Endo Optiks Service |
| All displays show “d” | When in pulse mode, the pulse exceeds the duration selected by 20%. | Recycle the AC power Call Endo Optiks Service |

Operator Replaceable Parts

The following list of spare parts is available from Endo Optiks or Endo Optiks distributors.

| <u>Part Name</u> | <u>Description</u> | <u>Part Number</u> |
|-------------------------|-----------------------------|--|
| Video Coupler | E2 MicroProbe Video Coupler | OME 300Z – 10k Scopes OME 300ZMG – 10k Scopes OME 170Z – 17k Scopes OME 170ZMG – 17k Scopes OME 600ZMG – 6k Scopes |
| Foot Switch | MicroProbe Foot Switch | SW001, SW002, SW003 |

TECHNICAL SPECIFICATIONS

Temperature

Celsius (C) and Fahrenheit (F)

- Operating temperature + 5° to 55° C (41° to 131 ° F)
- Storage (6 months) -30° to 50° C (-22° to 122° F)
- Transient (72 hours) -40° to 65° C (-40° to 149° F)

Humidity

- Operating 5% to 90% noncondensing
- Storage 5% to 95% noncondensing

Altitude

- Meters (m) and feet (ft.)
- Operating -305 to 3,048 m
(-1000 to 10,000 ft.)
- Shipping, non-operating -305 to 15,240 m
(-1000 to 50,000 ft)

Voltage requirements

- Voltage 115 AC or 240 AC

Power consumption (typical)

- Watts 690W or 720W

TECHNICAL SPECIFICATIONS

Light Source

175 Watt Xenon Light

- Visible Output 2,200 Lumens
- Radiant Output 25 Watts

300 Watt Xenon Light

- Visible Output 5,000 Lumens
- Radiant Output 50 Watts

Expected lifetime of the light source is 1,000 hours until the light reaches half of its initial output. (This could be equivalent to up to 6 years of usage.

Calculations are based on 3.25 hours of use per week.)

The failure of the light source is expected to be the deciding factor of the service life of the equipment.

CAUTION: The heat generated by the light source varies and can be hot causing the light connector, on the endoscope, to become hot. To prevent injury to the operator the light intensity should be turned down so the connector can cool before removing it from the system.

Treatment Laser

- Semiconductor GaAlAs Laser Diode, 810 ± 25 Nanometers,
- Waveform can be continuous or intermittent
- 1.2 Watts Output
- Class IV Laser
- Nominal Ocular Hazard Distance - 0.6 Meters

Aiming Beam Laser

- Visible Laser Diode, 640 ± 30 Nanometers,
- 1.5 Milliwatt Output
- Class IIIR Laser

Operational Mode

- Continuous Wave

Laser Source Power

- Less than 5 watts

Endoscope

- Beam Divergence - 40°
- Laser Delivery System - 200 μm fiberoptic terminating in a 0.89 mm diameter distal end SS tube with a laser fiber, video fiber and light fiber inside.
- Fiberoptic Connector - SMA 905 with the Endo Optiks Optical Safety Interlock to prevent firing the laser without the Endoscope connected.
- Nominal Ocular Hazard Distance - 0.6 Meters

Console Cooling

- Forced air

EMC TECHNICAL INFORMATION

The E4 Endoscopy System needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following information.

| Guidance and manufacturer's declaration - electromagnetic emissions | | |
|--|-------------------|---|
| <p>The E2 Microprobe Laser and Endoscopy System is intended for use in the electromagnetic environment specified below. The customer or the user of the E2 Microprobe Laser and Endoscopy System should assure that it is used in such an environment.</p> | | |
| Emission test | Compliance | Electromagnetic environment - guidance |
| RF emissions CISPR 11 | Group 1 | The E2 Microprobe Laser and Endoscopy System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class A | <p>The E2 Microprobe Laser and Endoscopy System is suitable for use in all establishments other than domestic ,and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment is intended for use by healthcare professionals only. This equipment may cause radio interference or may disrupt the operation of nearby equipment. It mat be necessary to take mitigation measures, such as re-orienting or relocating the E2 Microprobe Laser and Endoscopy System or shielding the location.</p> |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Complies | |

Guidance and manufacturer's declaration - electromagnetic immunity


The E2 Microprobe Laser and Endoscopy System is intended for use in the electromagnetic environment specified below. The customer or the user of the E2 Microprobe Laser and Endoscopy System should assure that it is used in such an environment.

| IMMUNITY test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|---|--|-------------------------|--|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6kV contact ±8kV air | compliant | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst | ±2kV for power supply lines ±1kV line(s) to earth | compliant | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ±1kV line(s) to line(s) ±2kV Line(s) to earth | compliant | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines. IEC 61000-4-8 | <5% <i>Ut</i> (>95% dip in <i>Ut</i>) for 0.5 cycle 40% <i>Ut</i> (60% dip in <i>Ut</i>) for 5 cycles 70% <i>Ut</i> (30% dip in <i>Ut</i>) for 25 cycles <5% <i>Ut</i> (>95% dip in <i>Ut</i>) for 5s | compliant | Mains power quality should be that of a typical commercial or hospital environment. If the user of the E2 Microprobe Laser and Endoscopy System requires continued operation during power mains interruptions, it is recommended that the E2 Microprobe Laser and Endoscopy System be powered from an uninterruptible power supply or battery. |
| Power frequency (50/60Hz) magnetic field IEC 61000-4-8 | 3A/m | compliant | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |

NOTE: *Ut* is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity

The E2 Microprobe Laser and Endoscopy System is intended for use in the electromagnetic environment specified below. The customer or the user of the E2 Microprobe Laser and Endoscopy System should assure that it is used in such an environment.

| IMMUNITY TEST | IEC 60601 TEST LEVEL | COMPLIANCE LEVEL | ELECTROMAGNETIC ENVIRONMENT - GUIDANCE |
|--|---|---------------------------|---|
| <p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p> | <p>3 Vrms 150kHz to 80MHz</p> <p>3 V/m 80MHz to 2.5GHz</p> | <p>[3]V</p> <p>[3]V/m</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the E2 Microprobe Laser and Endoscopy System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = [3.5/3]\sqrt{P}$</p> <p>$d = [3.5/3]\sqrt{P}$ 80MHz to 800MHz $d = [7/3]\sqrt{P}$ 800MHz to 2.5GHz</p> <p>where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM, FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the E2 Microprobe Laser and Endoscopy system is used exceeds the applicable RF compliance level above, the E2 Microprobe Laser and Endoscopy System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the E2 Microprobe Laser and Endoscopy System.
- b Over the frequency range 150kHz to 80MHz, field strengths should be less than [3]V/m.

Portable and mobile RF communications equipment can affect the E2 Microprobe Laser and Endoscopy System.

| Recommended separation distances between portable and mobile RF communications equipment and the E2 Microprobe Laser and Endoscopy System | | | | |
|--|--|---|--|----------|
| The E2 Microprobe Laser and Endoscopy System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the E2 Microprobe Laser and Endoscopy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the E2 Microprobe Laser and Endoscopy System as recommended below, according to the maximum output power of the communications equipment. | | | | |
| Rated maximum output power of transmitter W | Separation distance according to frequency of transmitter | | | m |
| | 150kHz to 80MHz $d = [1.667]\sqrt{P}$ | 80MHz to 800MHz $d = [1.667]\sqrt{P}$ | 800MHz to 2.5GHz $d = [2.333]\sqrt{P}$ | |
| 0.01 | d = 0.167 | d = 0.167 | d = 0.233 | |
| 0.1 | d = 0.527 | d = 0.527 | d = 0.738 | |
| 1 | d = 1.667 | d = 1.667 | d = 2.333 | |
| 10 | d = 5.271 | d = 5.271 | d = 7.377 | |
| 100 | d = 16.67 | d = 16.67 | d = 23.333 | |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. | | | | |
| NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies. | | | | |
| NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from the structures, objects and people. | | | | |

BIBLIOGRAPHY- ENDOSCOPY

1. Thorpe HE. Ocular endoscope. *Trans Am Acad Ophthalmol, Otolaryngol* 1934;39:422-4.
2. Norris JL, Cleasby GW. An endoscope for ophthalmology. *Am J Ophthalmol* 1978;85:420-2227.
3. Norris JL, Cleasby GW. Intraocular endoscopic surgery. *Am J Ophthalmol* 1981; 91:603-6.
4. Norris JL. Vitreous surgery viewed through an endoscope. *Dev. Ophthalmol.* 1981;2:15-16.
5. Norris JL, Cleasby GW. Intraocular foreign body removal by endoscopy. *Ann Ophthalmol.* 1982;14:371-372.
6. Shields MB. Intraocular cyclophotocoagulation. Histopathologic evaluation in primates. *Arch Ophthalmol* 1985;103:1731-3529.
7. Shields, M.B. Intraocular Cyclophotocoagulation *Trans.Ophthalmology Society of the U.K.* 1986; 105: 237-241.
8. Lecoq PJ, Billotte C, Combe JC. Interet de la videoendoscopie vitreo-retinienne. *J Fr Ophthalmol.* 1986;9:427-429.
9. Furia M, Hamard H, Puech M. Endoscopie oculaire. I. Modele experimental d'etude de l'implantation en chambre posterieure apres extraction extra-capsulaire du cristallin. *Bull. Soc. Ophthalmol Fr.* 1987;87:759-760.
10. Furia M, Hamard H, Puech M, Despreaux C, Elalouf M. Cloutage Retinien avec film video-endoscopique U-MATIC. *Bull Soc. Ophtalmol Fr.* 1987;87:1395-1403.
11. Lecoq PJ, Billotte C, Combe JC, Hamel C. Plaidoyer en faveur de l'endoscopie pour certaines interventions retinovitreennes. *Bull Oc Ophtalmol Fr.* 1987;87:575-576.
12. Kora Y, Yaguchi S. Sutured secondary posterior chamber lens with endoscopic control. *Ocular Surgery News,* 1990.
13. Volkov VV, Danilov AV. Flexible endoscope for intraocular surgery. *Arch Ophthalmol* 1990;108:1037-8
14. Volkov VV, Danilov AV. Flexible endoscopes: Ophthalmoendoscopic techniques and case reports. *Arch Ophthalmol* 1990;108:956-7
15. Eguchi S, Araie M. A new ophthalmic electronic video endoscope system for intraocular surgery. *Arch Ophthalmol* 1990;108:1778-81
16. Leon CS, Leon JA. Endoscopie Chirurgicale Oculaire. Medsi/McGraw-Hill, Paris France: 1990.
17. Leon CS, Leon JA. Microendoscopic ocular surgery Part II. Preliminary results from the study of glaucomatous eyes. *Journal of Cataract and Refractive Surgery* 1991;5:573-6
18. Kuhn F, Witherspoon CD, Morris RE. Endoscopic surgery vs temporary keratoprosthesis vitrectomy. *Arch Ophthalmology* 1991 Jun;109(6):768.
19. Uram M. Ophthalmic laser microendoscopy ciliary process ablation in the management of neovascular glaucoma. *Ophthal* 1992;99:1823-8.
20. Uram M. Ophthalmic laser microendoscope endophotocoagulation. *Ophthalmol* 1992;99:1829-32
21. Levin PS, Stormogipson J. Endocanalicular laser-assisted dacryocystorhinostomy. *Arch Ophthal* 1992;110:1488-90
22. Singh AD, Singh A, Whitmore I, Taylor E. Endoscopic visualization of the human nasolacrimal system: an experimental study. *BJO* 1992;76:663-7
23. Uram, M. Transcleral Nd:YAG laser cyclodestruction. *Ophthalmic Surgery* 1993;24:133
24. Joos MJ, Alward WLM, Folberg R. Experimental endoscopic goniotomy. *Ophthalmol* 1993;100:1066-70
25. Althaus C, Sundmacher R. Endoscopically controlled optimization of trans scleral suture fixation of posterior chamber lenses in the ciliary sulcus. *Ophthalmologie* 1993; 90:317-24
26. Joos KM, Shen J-H, Parel J-M, Rol P. In vitro examination of the anterior chamber angle with a gradient-index (GRIN) lens endoscope. *Lasers in Ophthalmology.* 1994;II 2330.
27. Fisher YL, Slakter JS. A Disposable Ophthalmic Endoscopic System. *Arch Ophthalmol.* 1994;112:984-986.

28. Kong YT, Kim TI, Kong BW. A report of 131 cases of endoscopic laser lacrimal surgery. *Ophthalmology* 1994; 101:1793-1800
29. Uram M. Laser endoscope in the management of proliferative vitreoretinopathy. *Ophthalmology* 1994; 101:1404-1408
30. Uram M. Diode laser endocycloablation. *Ophthalmic Surg* 1994;25:268-269
31. Mizota A, Takaso M, Asangi K, et al. Internal contact sclerostomy with an erbium laser and intraocular fiberscope. *Laser light Ophthalmol* 1995; 7:57 – 64
32. Trevisani MG, Allingham RR, Shields MB. Histologic comparison of contact transcleral diode cyclophotocoagulation and endoscopic diode cyclophotocoagulation. *Investigative Ophthalmol and Visual Science* 1995; 36:4
33. Uram M. Endoscopic cyclophotocoagulation in glaucoma management. *Current Opinion in Ophthalmology* 1995; 11:19-29
34. Uram M. Endoscopic Cyclophotocoagulation in Glaucoma Management: Indications, Results, and Complications. *Ophthalmic Practice* 1995; 13:173-185
35. Uram M. Combined Phacoemulsification, Endoscopic Cyclophotocoagulation, and Intraocular Lens Insertion in Glaucoma Management. *Ophthalmic Surgery* 1995; 26:346-352
36. Uram M. Transcanalicular diode laser dacryocystorhinostomy. *American Society of Cataract and Refractive Surgery*, June 1-4, 1996, Seattle, WA.
37. Gayton, J. Combined surgery with an endoscopic laser. *Maximizing Results: Strategies in refractive, corneal, cataract and glaucoma surgery*. Pp.224-234 Thorofare, N.J. SLACK Inc; 1996
38. Uram M. Endoscopic fluorescein angiography of the ciliary body in glaucoma management. *Ophthalmic surgery and lasers* 1996; 27:174-178
39. Uram M. Endoscopic Fluorescein Angiography. *Ophthalmic surgery and lasers* 1996;27:849-855.
40. Uram M. "Endoscopic Cyclophotocoagulation in the management of Aphakic and Pseudophakic Glaucoma" *Investigative Ophthalmol and Vis Sci* 1996;37:262
41. Jacobi PC, Dietlein TS, Krieglstein GK. Experimental microendoscopic photoablative laser goniotomy as a surgical model for the treatment of dysgenetic glaucoma. *Graefes Arch Clin Exp Ophthalmology*. 1996 Nov;234(11):670-6.
42. Jurgens I, Lillo J, Buil JA, et.al. Endoscope assisted transcleral suture fixation of intraocular lenses. *J Cataract Refract Surg* 1996;22:879-81
43. Tsunoda K. Treatment of anterior vitreous before suturing an intraocular lens to the ciliary sulcus. *J Cataract Refract Surg* 1996; 22:222-6
44. Terasaki H, Miyake Y, Awaya S. Fluorescein angiography of peripheral retina and pars plana during vitrectomy for proliferative diabetic retinopathy. *Am J Ophthalmology*. 1997 Mar;123(3):370-6
45. Medow N, Haley J, Lima F. Initial ciliary ablation with TSPC [letter]. *Ophthalmology* 1997;104:171-172
46. Mora JS, Iwach AG, Gafney MM, et. Al. Endoscopic diode laser cyclophoto-coagulation with a limbal approach. *Ophthalmic surg lasers* 1997; 28:118-123
47. Chen J, Cohn RA, Lin SC, et al. Endoscopic photocoagulation of the ciliary body for treatment of refractory glaucomas. *Ophthalmol* 1997; 124:787-796
48. Koch FHJ, Luloh KP, Augustin AJ, et. al. Subretinal surgery with gradient index endoscopes. *Ophthalmologica* 1997; 211:283-287
49. Leagis JM, Rol P, Briat B, et. al. Endoscope rigide a lentilles de GRIN. *J Fr Ophthalmology* 1997; 20:439-443
50. Medow NB, Sauer HL. Goniotomy for Congenital Glaucoma. *Journal Ped. Ophth. & Strabismus* 1997; 34:258-259.
51. Boscher C., Lebuissou DA, Jean JS, Nguyen-Khoa JL. Vitrectomy with endoscopy for management of retained lens fragments and/or posteriorly dislocated intraocular lens. *Graefes Arch Clin Exp Ophthalmol*. 1998 Feb;236(2):115-21.

52. Wallace DK, Plager DA, et al. Surgical Results of Secondary Glaucomas in Childhood. *Ophthalmology* 1998;105:101-111.
53. Scheindlin JA, Hirose T. Ophthalmic Endoscopy. *IOVS* 1998; 39:480.
54. Caronia RM, Sturm RT, Marmor MA, et. al. Treatment of a cyclodialysis cleft by means of ophthalmic laser endoscope endophotocoagulation. *AJO* 1999; 760-61
55. Gayton JL. Combined cataract and glaucoma surgery: trabeculectomy vs. endoscopic laser cyclophotocoagulation. *J ASCRS* 1999; 25:1212-1219
56. Feilder AR, Bentlet C, Moseley MJ. Recent Advances in Ophthalmology. *BMJ* 1999; 318:717-720.
57. Jacobi PC, Dietlein TS, et. al. Microendoscopic Trabecular Surgery in Glaucoma Management. *Ophthalmology* 1999; 106:538-544.
58. Scheindlin JA, Hirose T, Hartnett ME. Ophthalmic endoscopy: applications in intraocular surgery. *Int Ophthalmol Clin.* 1999 Winter;39(1):237-47.
59. Koch FH, Gumbel HO, Hattenbach LO, Ohrloff C. Intravitreal endoscopic visualization of intraocular ganciclovir devices: improved long-term treatment of CMV retinitis. *Klin Monatsbl Augenheilkd.* 1999 Feb;214(2):107-11
60. Plager DA, Neeley DE. Intermediate-term results of endoscopic diode laser cyclophotocoagulation for pediatric glaucoma. *J AAPOS* 1999; 3:131- 137
61. Uy HS, Foster CS. Diagnostic vitrectomy and uveitis. *Int. Ophthalmol Clin* 1999; 39:223-235.
62. Jacobi PC, Dietlein TS. Endoscopic surgery in glaucoma management. *Current Opinion in Ophthalmology.* 2000 Apr;11(2):127-32.
63. Marc D. DE SMET. Retinal Transscleral Photocoagulation under endoscopic control. *Retina* 2000; 20:3:316-317
64. Berke SJ, Cohen AJ, Sturm RT, Nelson, DB. Endoscopic cyclophotocoagulation and phacoemulsification in the treatment of medically controlled primary open-angle glaucoma. *Journal of Glaucoma* 2000;9:129.
65. Jacobi P. Deitlin T. Endoscopic surgery in glaucoma management. *Curr Opinion in Ophthalmol* 2000; 11: 127-132
66. Joos KM, Shen JH. An Ocular Endoscope Enables a Goniotomy Despite a Cloudy Cornea. *Arch. Ophthalmology* 2001;119:134-135.
67. Ciardella AP, Fisher YL, et. al. Endoscopic Vitreoretinal Surgery For Complicated Proliferative Diabetic Retinopathy. *Retina* 1/2001;21:20-27.
68. Singh AD et al. Endoscopic visualization of the human nasolacrimal system: an experimental study. *Brit. Journal of Ophthalmology* 1992;76:663-667.
69. Ben-nun J., Cornea sparing by endoscopically guided vitreoretinal surgery. *Ophthalmology* 2001;108:1465-1470.
70. Neely, Plager. Endocyclophotocoagulation for management of difficult pediatric glaucomas. *Journal of AAPOS* 2001;5:221-229.
71. Tello C, Walsh J. Rosen R, et. al. High resolution endoscopic evaluation of anterior segment structures. *Investigative ophthalmology and visual science* 2001; 42:767.
72. Agarwal A., Brummer TA, Tawansy KA, et al. Endoscopic assisted pars plana vitrectomy and glaucoma shunt placement. *Investigative ophthalmology and visual science* 2001;42:2349.
73. Torii H. Yoshitomi K, Miyata K., et al. Mechanical detachment of the anterior hyaloid membrane from the posterior lens capsule. *Ophthalmology* 2001;108:2182-2185.
74. Barraktar S, Koseoglu T. Endoscopic goniotomy with anterior chamber maintainer: Surgical technique and one year results. *Ophthalmic surg and lasers* 2001: 32;496-502
75. Kawai K. The microendoscope for ciliary process photocoagulation in neovascular glaucoma. *Tokai J Exp Clin Med.* 2002 Apr;27(1):27-34.
76. Barkana Y, Morad Y, Ben-nun J. Endoscopic photocoagulation of the ciliary body after repeated failure of transcleral diode laser cyclophotocoagulation. *American Journal of Ophthalmology* 2002;133:405-407.
77. Lin S. Endoscopic cyclophotocoagulation. *British Journal of ophthalmology* 2002;86:1434-38.

78. Sheilds S., Chen P. Sequential or simultaneous cyclophotocoagulation and glaucoma drainage implant for refractory glaucoma. *J Glaucoma* 2002; 11:203-208
79. Mizota A, Takasoh M, Kobyashi K, et. al. Internal sclerostomy with the Er:YAG laser using a gradient index lens endoscope. *Ophthalmic Surg Lasers* 2002; 33:214-20
80. Valtot F. The refractory glaucomas. *J Fr Ophthalmol.* 2003 Oct;26 Spec No 2:S56-61.
81. Huang S, Yu M. Application of endoscopic operation in ophthalmology. *Yan Ke Xue Bao.* 2003 Dec;19(4):227-31,243.
82. Hammer ME, Grizzard S. Endoscopy for evaluation and treatment of the ciliary body in Hypotony. *Retina* 2003;23:33-36.
83. Uram M. Endoscopic Surgery in Ophthalmology. 2003. Lippincott Williams Wilkens, Philadelphia, PA.
84. Rhee D. Endoscopic cyclophotocoagulation. Color Atlas and Synopsis of Clinical Ophthalmology: Glaucoma. 2003. McGraw-Hill, New York, New York
85. Hollander DA, Lin SC. Delayed therapeutic success with endoscopic cyclophotocoagulation. *BJO* 2003; 87:792-793
86. Faude F., Wiedemann P. Vitreoretinal endoscope for the assessment of the peripheral retina and the ciliary body after large retinectomies in sever anterior retina and the ciliary body after large retinectomies in severe anterior PVR. *Int. Ophthalmol.* 2004 Jan;25(1):53-6.
87. Schlote T. Surgical management of posttraumatic secondary angle closure glaucoma. *Klin Monatsbl Augenheilkd.* 2004 Aug;221(8):646-51.
88. Gierek-Lapinska A, Leszczynski R. Laser therapy in the treatment of glaucoma. *Klin Oczna.* 2004;106(1-2 Suppl):269-72.
89. Rosen KD, Chynn EW. Endoscopic surgery in ophthalmology. Book review. *Ocular Surgery News* 1/1/2004.
90. Valmaggia C, de Smet M. Endoscopic laser coagulation of the ciliary processes in patients with severe chronic glaucoma. *Klin Monatsbl Augenheilkd.* 2004 May; 221(5):343-6.
91. Moore JE, Herath GD, Sharma A. Continuous curvilinear capsulorhexis with use of an endoscope. *Journal of Cataract & Refractive Surgery.* May 2004; 30:960-963.
92. Koch FH, Quiroz-Mercado H, Hattenbach LO, Guerra-Naranjo JL, Yesuran I, Augustin A, Luloh KP, Ohrloff C. Pigment epithelium endoscopic laser surgery for treatment of choroidal neovascularization. *Ophthalmology.* 2004 May-June;218(3):162-75.
93. Lima FE, Magacho L, Carvalho DM, Susanna R Jr, Avila MP. A prospective, comparative study between endoscopic cyclophotocoagulation and the Ahmed drainage implant in refractory glaucoma. *J Glaucoma* 2004 Jun;13 (3):233-7.
94. Lin M, Huang S, Yu M. The Histopathologic changes of human eyes after laser endoscope cyclophotocoagulation. *Yan Ke Xue Bao* 2004 Dec; 20:233-6
95. Faude F, Weidemann P. Vitreoretinal endoscope for the assessment of the peripheral retina and ciliary body after large retinectomies in severe anterior PVR. *International Ophthalmology* 2004; 25:53-56
96. Fu AD, McDonald HR, Jumper MJ, et al. Recurrent vitreous hemorrhage after sutured posterior chamber intraocular lenses. *Retina* 2004; 24:193-198
97. Poothulli AM, McLeod SD, Lin S. Employing endoscopic guidance for placement of a black diaphragm aniridia intraocular lens following destructive *Acanthamoeba* sclerokeratitis. *British J Ophthalmol* 2004; 88:596-597
98. Holz HA, Lim MC. Glaucoma lasers: a review of newer techniques. *Curr Opin Ophthalmol* 2005 April; 16:89-93
99. Uka J, Minamoto A, Hirayama T, et. al. Endoscope aided cataract surgery in corneal opacity associated with aniridia. *J Cataract and Refractive Surg* 2005; 31:1455-57
100. Hong JE, Hatton HP, Leib M, et. al. Endocanalicular laser dacryocystorhinostomy. *Ophthalmology* 2005;112:1629-1633.

101. Chen SDM, Salmon JF, Patel CK. Video endoscope guided fluorescein assisted vitrectomy for phakic malignant glaucoma. *Arch Ophthalmol* 2005;123:1419-1421
102. Carter CB, Plager DA, Neely DE, et. al. Endoscopic diode laser cyclophotocoagulation in the management of aphakic and pseudophakic glaucoma. *JAAPOS* 2006;10:63-73
103. Kleinmann G, Apple DJ, Chew J, Mamalis N. New Endoscopic Technique to Analyze Various Modern Specialized Intraocular Lenses in Research Eyes and Human Eyes Obtained Postmortem. *Ophthalmology* April 2006; 113 (4): 591-596.
104. Berke SJ. Endolaser Cyclophotocoagulation in Glaucoma Management. *Techniques in Ophthalmology* 4(2):74-81, 2006.
105. Bloom PA, Dharmaraj S. Endoscopic and transscleral cyclophotocoagulation. *BJO* 2006;000:1-2.
106. Lin SC, Chen MJ, Lin MS, et al. Vascular effects on ciliary tissue from endoscopic versus transscleral cyclophotocoagulation. *Br J Ophthalmol* 2006;90:526-30.
107. Sonoda Y, Yamakiri K. Endoscopy guided subretinal fluid drainage in vitrectomy for retinal detachment. *Ophthalmologica* 2006;220:83-86
108. Henderson B, Ali R, Kim JY, Shortsleeve Ament C. Using endoscopy to teach cataract surgery. *J Cataract Refract Surg.* 2006;032: 1606-1610.
109. Davis NL, Wetli CV, Shakin JL. The retina in forensic medicine. Applications of ophthalmic endoscopy: the first 100 cases. *Am J for Med Path* 2006:27-10.
110. Yu MB, Huang SS, Ge J, et. al . The clinical study of endoscopic cyclophotocoagulation on the management of refractory glaucoma. *Zhonghua Yan Ke Za Zhi* 2006; 42:27-31
111. Pantcheva MB, Kahook MY, Schuman JS, Noecker RJ. Comparison of acute structural and histopathologic changes in human autopsy eyes after endoscopic cyclophotocoagulation and tran-scleral cyclophotocoagulation. *Br J Ophthalmol* 2007; 91 : 248-252
112. Kahook MY, Schuman JS, Noecker RJ. Endoscopic cyclophotocoagulation using iris hooks versus viscoelastic devices. *Ophthalmic Surg lasers imaging* 2007;38:170-172
113. Pantcheva MB, Kahook M, Schuman JS, et al. Comparison of acute structural and histological changes of the porcine ciliary processes after endoscopic cyclophotocoagulation and transcleral cyclophotocoagulation. *Clin Experiment Ophthalmol.* 2007; 35:270-274
114. Kahook MY, Lathrop KL, Noecker RJ. One site versus two site endoscopic cyclophotocoagulation. *Journal of Glaucoma* 2007;16:527-530 MY, Lathrop KL, Noecker RJ. One site versus two site endoscopic cyclophotocoagulation. *Journal of Glaucoma* 2007;16:527-530
115. Peyman G. *Vitreoretinal Surgical Techniques.* 2007. Taylor and Francis, New York, NY
116. Ramula PY, Corcoran SL, et. al. Utilization of Various Glaucoma Surgeries and Procedures in Medicare Beneficiaries from 1995 to 2004. *Ophthalmology* 2007;114:2265-2270
117. Minckler DS, Francis BA, Hodapp EA, et. al. Aqueous shunts in glaucoma. *Ophthalmology* 2008;115:1089-1098 (TUBES AND TRABS FAIL @ 10%/YR)
118. Parthasarathy A, Aung T, Oen FT, Tan DT. Endoscopic cyclophotocoagulation for the management of advanced glaucoma after osteo-odonto-keratoprosthesis surgery. *Clin Experiment Ophthalmol* 2008;36:93-94
119. Lee SM, Kim mk, Oh JY, et. al. Endoscopic vitrectomy improves outcomes of Seoul-type keratoprosthesis exchange in rabbit model. *Investigative Ophthalmol and Vis Sci* 2008;49:4407-4411
120. Weston J. Techniques of Endoscopic Cyclophotocoagulation. *Techniques in Ophthalmology* 6(3):98-104, 2008
121. Al-Haddad CE, Freeman SF. Endoscopic laser cyclophotocoagulation in pediatric glaucoma with corneal opacities. *JAAPOS* 2008;11:23-8
122. Nagpal M, Wartikar S, Nagpal K. Endoscopic evaluation of sclerotomy sites. *Retina* 2009;29:225-231
123. Murthy GJ, Murthy PR, Murthy KR, et. al. A study of the efficacy of endoscopic cyclophotocoagulation for the treatment of refractory glaucomas. *Indian J Ophthalmol.* 2009;57:127-132

124. Falkenberry SM, Siegfried C. Endocyclophotocoagulation. *Middle East Afr Ophthalmol.* 2009;16:130-133
125. Vold SD. , Cataract surgery combined with endoscopic cyclophotocoagulation. *Cataract Surgery in the Glaucoma Patient* , Springer New York, 2009, pp129-134
126. Henderson B, Ali R, Kim J., et al. Using endoscopy to teach cataract surgery. *J Cataract and Refr Surg* 2009; 32:1606-1610
127. Fine IH, Mojon DS. Minimally Invasive Ophthalmic Surgery. Springer Berlin Heidelberg; 2010, p.217-231
128. Fang A, Yang X, Nie L, et. al. Endoscopically controlled goniosynerchysis in managing synechial angle closure glaucoma. *J Glaucoma* 2010; 19:19-23
129. Kawashima M, Kawashima S, Dogru M, et.al. Endoscopy guided vitreoretinal surgery following a penetrating corneal injury. *Clinical Ophthalmology* 2010; 4:895-898
130. Haverly R. Surgical management of neovascular glaucoma. *Tech Ophthalmol* 2010; 8:87-93.
131. Francis, Kawji, Vo, et.al. Endoscopic Cyclophotocoagulation (ECP) in the Management of Uncontrolled Glaucoma with prior aqueous tube shunt. *J. Glaucoma* 2010; 00-000-000
132. Raizada S, Sabti K. Endoscope assisted enhanced internal search for iatrogenic retinal breaks in 20 gauge macular surgery. *Br J Ophthalmol* 2010; 10:1136.
133. Olsen TW, Pribila, JT. Pars Plana Vitrectomy With Endoscope-Guided Sutured Posterior Chamber Intraocular Lens Implantation in Children and Adults. *AJO* 2011; 151:287-296
134. Tarantola, RM., Agarwal, A., et.al. Long-term results of combined endoscope-assisted pars plana vitrectomy and glaucoma tube shunt surgery. *Retina* 31:275-283, 2011.
135. Kita, M., Yoshimura, N. Endoscope-Assisted Vitrectomy in the Management of Pseudophakic and Aphakic Retinal Detachments with Undetected Retinal Breaks. *Retina* 0:1-5, 2011.

**HERE ARE SOME SUGGESTIONS
REGARDING ECP CASE SELECTION AND
EXTENT OF TREATMENT!**

| | | |
|--|-------------|---|
| COMBINED WITH PHACO WITH MEDICALLY CONTROLLED GLAUCOMA PRE-OPERATIVELY | 180° | GOOD FOR LEARNING |
| WITH UNCONTROLLED GLAUCOMA | 270° - 360° | GOOD FOR LEARNING |
| POSTERIOR CHAMBER PSEUDOPHAKIA | | |
| OPEN-ANGLE GLAUCOMA | 270° - 360° | GOOD FOR LEARNING |
| NON-VASCULARIZED ANGLE CLOSURE | 270° - 360° | GOOD FOR LEARNING |
| ANTERIOR CHAMBER PSEUDOPHAKIA | 270° - 360° | CAN BE DIFFICULT |
| APHAKIC OPEN OR CLOSED ANGLE GLAUCOMA | 270° - 360° | CAN BE DIFFICULT |
| PHAKIC OPEN OR CLOSED ANGLE GLAUCOMA | 270° - 360° | SHOULD HAVE SOME PRIOR EXPERIENCE |
| NEOVASCULAR GLAUCOMA | 180° FIRST | MAY NEED ADDITIONAL 90° - 180° LENS STATUS DETERMINES DIFFICULTY |

QUICK SET UP GUIDE

Color Video Monitor

⇒ Turn monitor ON by depressing power button **lower right**

Back of E2 Unit

⇒ Start with Laser Key in the OFF Position (O) **top center**

⇒ Turn Power ON (1) **above power cord**

⇒ Laser Key to ON Position (1)

E2 Front Control Panel

⇒ Connect Three Probe Ends to Front Panel (save caps)

Laser - Screw on - **top right**

Light - Plug in - **lower center**

Video Connection - Firm Plug in - **lower right**

⇒ Laser Power (**yellow**) - set to .25 watts

⇒ Laser Duration (**blue**) - set to **Continuous**

⇒ Increase Illumination (**white**) - only partially

⇒ Focus Image - knurled ring around video connection.

⇒ Press Enable Button (**green**) - indicator light comes on

⇒ Increase Aiming Beam (**red**) - red light in field

- When the surgeon is ready to use microendoscope, light and aiming beam may need to be adjusted.
- Laser power may be adjusted intraoperatively while pedal is not depressed.
- Red emergency stop button must be in the **out** (off) position for laser to operate.
- ONLY sodium hyaluronate viscoelastics may be instilled when laser is to be applied.
- If image becomes blurred during surgery, clean probe tip with Meroceel® Instrument Wipe or Spear.
- Before disconnecting microendoscope, clean probe tip with concentrated alcohol and then rinse. Check image clarity before turning off unit.

CAUTION: The heat generated by the light source varies and can be hot causing the light connector, on the endoscope, to become hot. To prevent injury to the operator the light intensity should be turned down so the connector can cool before removing it from the system.