


*Kenosha Medical Center  
SN-ES171M9811*

ENGLISH E  
OPERATOR GUIDE  
*Cont. #*  
**77210163 GB**

CRYOSTAT SERIES





ENGLISH E

# OPERATOR GUIDE

CGH = #  
**77210163 GB**



Life Sciences International (Europe) Limited is an ISO 9001 Approved Company.

# WELCOME

## 1 INTRODUCTION

Welcome to the Cryotome® Electronic (E) cryostat - . Designed and made with care, the instrument is safe to use, simple to operate, and easy to maintain.

This Operator Guide gives instructions for its correct operation and use. The number by the side of each illustration is the same as the paragraph number of its associated text.

## 2 SAFETY

Shandon instruments are designed to operate safely and consistently for many years. However, incorrect actions by a user may damage the equipment, or cause a hazard to health. It is important for you to know that:

- i The instrument weighs approximately 150 Kilograms (330 lbs); get help to move or lift it.
- ii Sharp components are used in this instrument. Make sure you understand the correct methods for fixing and use.
- iii Be aware that there is always a danger of biological contamination when working with unfixed tissue. Make sure that the instrument is properly decontaminated at regular intervals.
- iv Do not remove any panels or covers while the instrument is attached to a mains supply.
- v This instrument must be properly connected to a good earth (Ground) via the mains input supply. All earth (ground) straps and connections **MUST** be replaced after service and before the Cryotome® is powered up.

In compliance with statutory requirements all our equipment is designed to accepted standards of safety. Its use does not entail any hazard if operated in accordance with the instructions given in this manual. However, you must obey the following safety precautions.

- i All users must read and understand the Operator Guide, and only operate the unit in accordance with the instructions.
- ii Potentially lethal voltages above 110v A.C. or 50v D.C. are present inside the unit. Do not remove any access covers unless specifically instructed to do so.
- iii It is important that normal standards of safety and good laboratory practices are employed, **especially with respect to decontamination and cleaning procedures**. Always use common sense when operating the instrument.
- iv Any problems and queries should be referred to our Service Department.
- v Do not use the Cryotome® for general refrigeration purposes, such as storage of specimens. Use a conventional freezer.
- vi Correct maintenance procedures are essential for consistent performance. It is recommended that a Maintenance Contract is taken out with our Service Department.
- vii It is important that only factory approved accessories or replacement parts are used with the Cryotome.

## 3 UNPACKING

Transit fixings are attached to various internal structures to protect the instrument during transport. **IT IS ESSENTIAL THAT THESE ARE ALL REMOVED BEFORE AN ATTEMPT IS MADE TO OPERATE THE INSTRUMENT.**

Refer to Section 2 for full instructions relating to unpacking and installation.

# DESCRIPTION

## 1.1 OVERVIEW

1.1.1 The Cryotome<sup>®</sup> is designed for cutting frozen sections of fixed and unfixed specimens, rapidly and accurately, for future examination by microscope. It operates by rapidly freezing the sample in a temperature controlled environment, and provides means for sectioning the frozen sample by conventional microtome. For optimum performance, the microtome Specimen Head is held within the refrigerated chamber. The main body of the microtome is mounted outside the chamber to give more space for working in, and for easier cleaning.

1.1.2 Principal factors for good sectioning of frozen specimens are:

- i the temperature must be correct for the specimen being cut.
- ii the microtome must be correctly adjusted and operated.
- iii the cutting blade must be sharp and set at the correct angle.
- iv the anti-roll plate must be correctly adjusted.


1.1.3 The Cryotome<sup>®</sup> meets these objectives by good design and high quality manufacture. In the Electronic version, additional sophistication is provided by electronic control of the specimen advance and retract functions.

1.1.4 The main body of the microtome outside the refrigerated chamber connects to the **SPECIMEN HEAD** via a thermal barrier that also serves to provide a biological screen.

TABLE 1

### RECOMMENDED TEMPERATURES FOR CUTTING UNFIXED FROZEN TISSUES.

TISSUE TYPE	WORKING TEMP.(°C)
Brain	-12
Liver	-14
Lymph Node	-14
Kidney	-16
Spleen	-16
Muscle	-20
Thyroid	-20
Skin	-25
Breast	-25
Breast with fat	-30 or below
Adipose tissue	-30 or below
Fixed Tissue	-12 to -17



1.1.4

1.1.9 Each type of Holder is specifically designed to accept one style of blade only. Both the Lo-Profile and Hi-Profile Holders are for use with disposable blades, and all knife-holders, feature an **ANTI ROLL PLATE** that prevents curling of the cut section.


1.1.10 The knife Holder is constructed of two parts - a Top and a Base. The Top holds the blade or knife, and the Base supports the Top. There are three types of Top to suit the types of knife or blade. The same Base supports all versions of Top.

1.1.11 The assembled Top and Base of the knife Holder fix onto the **CARRIAGE** on the floor of the Refrigerated Chamber. Four Allen screws hold the Carriage in place and attach it to the body of the microtome outside the Refrigerated Chamber. These Allen screws must be kept tightly fastened to prevent movement of the Holder - and hence the blade - with respect to the Specimen Head which is essential for good sectioning.


**WARNING** Do not undo or remove the Allen screws or the whole of the Microtome assembly will fall and cause considerable damage.

1.1.12 The **FLYWHEEL** on the righthand side of the instrument controls the movement of the Specimen Head, and enables the Cutting Sector to be set accurately. It also enables sections to be taken manually if required.


1.1.9




1.1.10



1.1.11



1.1.12




1.1.16 A trough, fitted behind the Cryobar on the rear wall of the Refrigerated Chamber accommodates approximately 2 ml of concentrated Formalin. Pressing **[IMMED FUMIGATE]** on the Control Panel, causes the trough to heat up and evaporate the Formalin during the course of a fumigation cycle. The formalin vapour (formaldehyde) decontaminates the working chamber. For safety, the fumigation cycle cannot take place if the window is not closed and locked.

1.1.17 The window is double glazed and contains an electrical de-misting element that is held permanently energised during normal operation of the Cryotome®. The window slides backward, under the Control Panel, to give access to the Refrigerated Chamber and is lockable when closed. A window-locked indicator on the Control Panel confirms this status.


1.1.18 A **WINDOW BOOST** switch next to the mains supply socket at the rear of the instrument provides facility for adding extra de-misting power in conditions of high humidity.

**WARNING** USE OF THE WINDOW BOOST FACILITY SHOULD BE KEPT TO A MINIMUM OR PERFORMANCE MAY BE ADVERSELY AFFECTED.

1.1.16



1.1.18



# INSTALLATION INSTRUCTIONS

## 2.1 PROCEDURE

2.1.1 Check that the detail on the label of the crate corresponds with the Purchase Order, and that the power supply capability of your local socket outlet is compatible with the power demand of the Cryotome.



**THE CRYOTOME IS HEAVY (150 Kgs 330 LBS) GET HELP TO MOVE OR LIFT IT.**

2.1.2 Move the crate near to where the instrument is to be sited. Cut the retaining straps then lift the outer case vertically. Read the label on the rear of the instrument and check that it conforms with your order.

2.1.3 Retrieve the accessories pack from on top of the pallet in front of the instrument, then lift the Cryotome off the Pallet. Check that the Accessories Pack contains:

- i Pack of Five Cryocassettes.
- ii One Shelf
- iii One Debris Tray
- iv One Flywheel Assembly
- v One Flywheel Bolt
- vi One Flywheel Bolt Key
- vii Mains Lead
- viii One Brush - for sections
- ix Two Keys - front Panel
- x Tools.
- xi Operator Guide.


## POWER DEMAND CHARACTERISTICS

NOMINAL A.C.INPUT	MINIMUM VOLTAGE @ I (max)	I (max) = AMPS WHEN COMPRESSOR STARTS	MAX VOLTS AT NO LOAD
230 50Hz	198V	20 Amps	264V
110 60Hz	99V	44 Amps	121
220 60Hz	198V	20 Amps	242

### Notes

- 1 I (max) = Peak Current
- 2 Peak Current is for 22 secs max with stalled rotor.
- 3 Voltage drop affects compressor start capability.
- 4 Supply impedance details are given on Page
- 5 50 Hz units = 1700 VA
- 6 60 Hz units = 1400 VA

### 2.1.2



2.1.7 Hold the blue plastic coil of the defrost valve firmly on the transit plate at the back of the microtome, then cut the ties that hold the valve .

2.1.8 With the coil of the defrost valve still connected to the white wires, place the coil over the tube that projects upward from the copper tube on top of the compressor cooling matrix. Release the blue plastic components and fit the nameplate, the washer, and then the nut to hold the coil of the defrost valve in position. Fit the pastic plug in the hole in the centre of the rear panel

2.1.9 From the underside, support the stud of the clamping bracket that prevents damage to the microtome leadscrew, then undo and remove the top nut. Support the stud from underneath while allowing it to fall clear of the leadscrew arm.


2.1.10 Use the nylon loop to remove the spacer from between the arm of the leadscrew and the casting of the microtome.

2.1.11 Check that the detection flag at the end of the arm is located between the two black plastic opto units of the printed circuit board. This signifies that the leadscrew and trunnion nut are located correctly.


**Note**

- 1 Turn the leadscrew to adjust the setting if necessary.


2.1.7




2.1.8




2.1.9



2.1.10




2.1.11




2.1.15 Raise the front edge of the top cover, remove the transit packing from the top of the sliding window, then lower the cover back into position. Push open the sliding window to gain access to the Refrigerated Chamber.

2.1.15




2.1.16 Remove the rubber band and any other packing material from the area around the Cryobar on the lefthand wall of the Refrigerated Chamber.

2.1.16



2.1.17 Slide the flywheel carefully onto the drive shaft at the top of the righthand side panel. The flywheel is supplied in the Accessory Pack. Make sure that the locating pin of the flywheel fits into the key-way of the drive shaft.

2.1.17




**Note**

- 1 *The lever next to the flywheel handle is used to lock the flywheel and prevent inadvertent rotation.*

2.1.18 Fit the flywheel bolt then use the special Flywheel Bolt Key to turn the retaining bolt clockwise to tighten. Rotate the flywheel and check that movement is smooth and produces a corresponding movement of the Specimen Head.

2.1.18





**MAKE SURE THAT BODY OF THE INSTRUMENT IS EFFECTIVELY EARTHED (GROUNDED). A PROPRIETARY GROUND FAULT CIRCUIT INTERRUPTION DEVICE (EARTH LEAK DETECTION) MAY BE FITTED AT THE WALL SOCKET AS AN ADDITIONAL SAFEGUARD**


2.1.23 Make sure that the **I / O** power switch at the rear of the instrument is set to off (**O** side pushed inward).

2.1.24 Connect the mains supply cable to the local power supply outlet. The Cryotome® is now fully installed, and ready to be set up for operation.


2.1.25 Press the **I (ON)** side of the **I / O** switch inward to switch the instrument on.

2.1.26 Insert the key in the keyswitch of the Control Panel and turn it clockwise. The **LOCKED** indicator is not lit when the panel is unlocked and available for use.

**2.1.23**



**2.1.26**



**WARNING** Practice good biological safety procedures whenever the Cryotome is used.

**WARNING** Make sure that the Waste Bottle is filled with sufficient 10% Formalin solution to cover the outlet pipe.

**WARNING** Do not use the Cryotome for general laboratory refrigeration purposes.

**WARNING** Fumigate regularly to ensure decontamination of the Refrigeration Chamber and accessories.

**WARNING** Use mesh gloves when handling knives and blades.

**WARNING** Use the correct type of knife or blade appropriate for the knife holder.

**WARNING** Make sure that the knife guard is correctly installed.

**WARNING** If you know that a specific virus or bacteria is likely to be present in the specimen, make sure that you are aware of a suitable decontaminant before you introduce the specimen into the Refrigerated Chamber.

**WARNING** Always make sure that the window is closed before any sections are cut. Aerosols can form when the specimen traverses the knife. These can be injurious to health if inhaled.

**THESE WARNINGS ARE INCLUDED FOR YOUR SAFETY. PLEASE MAKE SURE THAT YOU READ AND UNDERSTAND THEM BEFORE OPERATING THE INSTRUMENT.**



# CONTROLS

## 3.1 DESCRIPTION

3.1.1 The **MAIN CONTROL PANEL** comprises groups of touch sensitive wipe-clean membrane switches and associated liquid crystal displays (LCD). There is no 'feel' when a switch operates. Operation is confirmed by an audible tone, a change in the display, or by an integral pushbutton indicator.

3.1.2 After a push-button is pressed its indicator flashes yellow while the instrument is making adjustments. The indicator shows steady yellow when the required status is achieved. Display panels show instrument status and error messages.


3.1.3 The **MINOR CONTROL PANEL** on the left of the window has push-buttons that control the movement of the Specimen Head, as follows;

REWIND	∧
ADVANCE	∨
FAST REWIND	∧∧
FAST ADVANCE	∨∨


There is also a notched Electronic Rotary Control that advances or retracts the Specimen Head 8 microns per notch dependent on the direction of rotation - clockwise to advance; counter clockwise to retract.

3.1.4 The **MICRON SELECTOR** section of the Control Panel displays the thickness in microns of the sections being taken.

3.1.1



3.1.3



### 3.3 SPECIMEN TRAVEL

3.3.1 The vertical bar graph contains separate segments that light up sequentially as the Specimen Head advances from its maximum available travel position. The proportion of lit to unlit segments indicates the amount of travel still available at any time.

3.3.2 When the instrument is first powered up, the first segment flashes to indicate that the Specimen Head is in its maximum available travel (MAX) position.

3.3.3 As the Specimen Head advances at the end of each cutting stroke, the available travel gets used up such that when no more travel is available all the **SPECIMEN TRAVEL** segments are lit. An audible alarm operates when the end of travel is reached.

3.3.4 Pressing [**RESET**] lights the integral push-button indicator, returns the Specimen Head to its maximum available travel position, and causes the display to flash.

3.3.5 When the instrument is being reset, all other advance or rewind functions are overridden. This occurs each time the instrument is switched on, or [**RESET**] is pressed.


### 3.4 TIME CONTROLS

3.4.1 The time setting controls comprise:


- i three illuminating push-buttons.
- ii three non illuminating push-buttons,
- iii an LCD display panel.

3.4.2 The illuminating push-buttons flash while an adjustment is being made after being pressed. Real time shows in 24 hour clock format, and a flashing colon (:) shows that operation is normal.

#### 3.3.1



#### 3.4.1



DAY 0 = prevents the programmed function.

DAY 1-7 = defrosting starts at the time set on the day shown (set Monday as Day 1 to avoid confusion).

DAY 8 = defrosting starts at the set time **every day**

DAY 9 = STANDBY (PSAVE = power save)

3.6.2 Standby (PSAVE) allows the instrument to remain switched on in a defrosted state for up to 24 hours until a pre-determined time. It then returns to its normal operational state

3.6.3 **TO SET STANDBY (PSAVE) mode** press and hold [DEFROST] then select **DAY 9** Set the time that normal operation is required to resume. Release the push-buttons. To start (PSAVE) mode immediately press [IMMED DEFROST]. The instrument then runs its defrost cycle for 15 mins and goes to (PSAVE) on completion. You can use [CANCEL DEFROST/ FUMIGATE] to terminate the function before the cycle is complete. Fumigation can only be cancelled at specific times (Note 2, Page 3.7).

3.6.4 Press [CANCEL DEFROST/FUMIGATE] on the CHAMBER CONTROLS panel to resume normal operation after PSAVE .


**Note**

1 It can take 1 - 2 hours after (PSAVE) mode to cool the Refrigerated Chamber sufficiently for use. Make allowance for this when setting the resume time.


### 3.7 TO SET FUMIGATE TIMES

3.7.1 As for SET DEFROST (3.6.1 to 3.6.5) but read [FUMIGATE] for [DEFROST]. Press [IMMED FUMIGATE] instead of [IMMED DEFROST] to enter (PSAVE) mode. This starts the 5 hour FUMIGATE cycle before going into (PSAVE).

### 3.6.3



### 3.6.4



**Note**

- 1 The compressor will not restart until 3 minutes after it stopped.


3.8.5 If the DAY 9 option is selected when setting the defrost times, then the instrument enters STANDBY (PSAVE) mode when defrosting finishes.

3.8.6 Pressing [IMMED FUMIGATE] initiates the following sequence:


**Notes**

- 1 The window **MUST** be shut and locked before fumigation can start. If [IMMED FUMIGATE] is pressed when the window is not shut and locked, **Err 4** is displayed and the fumigate cycle is aborted.
- 2 [CANCEL DEFROST/FUMIGATE] only cancels FUMIGATE during stages iii and iv that follow. Stages v and vi cannot be cancelled.
  - i the [IMMED FUMIGATE] push-button indicator lights and the display shows **Ald** (Formaldehyde).
  - ii the window lock is activated if the window is not already locked by manual use of the [WINDOW LOCK] push-button.
  - iii a 5 minute defrost is started, then....
  - iv the compressor is switched off and a one hour defrost starts,
  - v the Formalin trough is heated to 120°C for 45 minutes
  - vi the heater is switched off and fumigation takes place over four hours.
  - vii the instrument reverts to its previous setting when the fumigate cycle is finished.

**3.8.6**



**3.8.6.i**



### 3.9 TEMPERATURE CONTROLS

3.9.1 Two illuminating push-buttons in the **TEMPERATURE °C** panel facilitate monitoring of component temperatures in the Refrigerated Chamber. They also operate in conjunction with non-illuminating push-buttons **[+]** and **[-]** to set temperatures.

3.9.2 Monitoring is effected by pressing and releasing **[CRYOBAR TEMP]** or **[WORKING TEMP]** as appropriate. The display then shows the actual temperature of that component in °C. A lit push-button shows the selection to which the display refers


3.9.3 To set the working temperature, press **[WORKING TEMP]** while using the **[+]** and **[-]** push-buttons to raise or lower the value of the number in the display. Press and hold until the desired temperature is displayed, then release both push-buttons.

#### Note


- 1 It is not possible to set the temperature of the Cryobar. The Cryobar is situated on the refrigeration pipework and its temperature can vary considerably.

3.9.4 The temperature in the Refrigerated Chamber (working temperature) is automatically controlled to within 1°C of the set temperature. If no temperature is set, or after a reset, the working temperature automatically defaults to -20°C, or the last programmed temperature.

3.9.5 **[CRYOBAR TEMP]** displays the surface temperature of the Cryobar. However, if **[CRYOBAR BOOST]** is lit on the adjacent panel, the display shows the approximate temperature of the small rectangular peltier element in the Cryobar.



3.9.3



3.9.5

- ii **Ald** shows when fumigation is in progress.
- iii **dEF** shows when defrosting is in progress.
- iv --- shows if the associated temperature probe is faulty.
- v **Fault Codes** show in the TEMPERATURE °C display if the instrument malfunctions. Press [**WORKING TEMP**] or [**CRYOBAR TEMP**] to over-ride most faults. Faults F:11, F:12 and F:13, cause the instrument to wait for 30 mins then reset and attempt to reach its default temperature of -20°C.

**Note**

1 If the fault fails to clear, contact your Shandon Dealer for service.

3.11.2 Indications that show in the **TIME CONTROLS** display are as follows:


- i **P FAIL** (flashing) shows after a normal switch OFF then ON, or if power failed and recovered. The TEMPERATURE °C panel shows the temperature of the Refrigerated Chamber at the time power was restored.

Press [**CLOCK**] to clear.


- ii **rESEt** (flashing) shows if a problem caused the instrument to reset and clear its memory to its default status. The number in TEMPERATURE °C is the temperature of the Refrigerated Chamber at the time normal operation resumed.

Press [**CLOCK**] to clear.


3.11.1(ii)




3.11.1(iii)



3.11.2 i



3.11.2 ii




# SETTING UP

## 4.1 INTRODUCTION


4.1.1 Setting up the Cryotome® for routine operation involves the installation and adjustment of one or more of the following components of the microtome.

- i Knife Holders
- ii Knife Guards.
- iv Standard solid Knife.
- v Cutting Angle.
- vi Disposable Blades.
- vii Disposable Blade cutting angle.
- viii Anti Roll Plate angle.
- ix Anti-Roll Plate parallelism.
- x Anti-Roll Plate gap.
- xi Anti Roll Plate height.

4.1.1(i)




4.1.1(ii)




4.2.3 Pull the righthand pivoting lever of the knife holder fully toward you.

4.2.3



4.2.4 Adjust the tilt of the knife holder until the '0' of the tilt scale is approximately in line with the marker on the knife base, then push the righthand lever away from you to lock.

4.2.4



**Note**

- i The knife holder is adjusted to its precise operating angle later, after the blade or knife is installed

### 4.3 TO FIT A STANDARD SOLID KNIFE




THE KNIFE IS SUBSTANTIAL, EXTREMELY SHARP, AND COULD CAUSE SERIOUS INJURY IF MIS-HANDLED. ALWAYS WEAR MESH GLOVES WHEN HANDLING A STANDARD SOLID KNIFE

**Note**

- 1 The knife guard is specifically designed to operate with 160X34X10mm knife profiles.

4.3.1 Loosen the two clamp knobs on the standard solid knife holder by turning them counter clockwise.

4.3.1



4.4.2 Swing the anti roll plate away from the blade clamp by turning the control bar toward you.


4.4.3 Move the blade clamp lever to the left (clockwise) until it is vertically down.

4.4.4 Push the end of the blade to be installed part way out of the dispenser and put its left end under the plate of the blade clamp. Slide the blade to the left until its right-hand end drops into the recess of the blade transporter.


4.4.5 Use the transporter to move the blade into its first cutting position, then push the blade clamp lever to the right (counter clockwise) to lock the blade in position.

4.4.6 With use it is possible that the blade clamp mechanism loosens slightly so that the blade does not clamp properly. Adjust the screw at the rear of the Knife Holder as necessary to obtain the required locking action.


4.4.2




4.4.3




4.4.4



4.4.5




4.4.6




4.5.5 Use the adjusting screw on the anti roll plate support bar to set the approximate height for the anti roll plate. Loosen and tighten the Allen screw as necessary.

4.5.5




4.5.6 Use the height adjust control to set the precise height of the anti roll plate. You should **JUST** be able to see the top of the anti-roll plate above the edge of the knife/blade. Fine adjust so that the specimen **JUST** misses the anti roll plate when cutting.

4.5.6



4.5.7 Use the nylon adjuster screw to rectify any slew of the anti roll plate. Make sure that the top of the anti roll plate is parallel with the top edge of the knife/blade.

4.5.7



4.5.8 Ensure that there is a good edge i.e. there are no 'nicks' - on the anti roll plate. Use the brush supplied to keep the anti roll plate, the blade, and the clamp clean.

4.6.4 To orientate the Specimen to cut the best sections, loosen the knob on top of the Specimen Head by turning it anti-clockwise. The Cryocassette Clamp may then rotate 360° in the cutting plane. Tighten the knob when the optimum angle is obtained.

4.6.5 Full orientation of the specimen is possible when a Fine Adjust Orientating Head (0620-006) is specified. This allows 360° rotation, as before, plus fine adjustment in the X-Y plane.

4.6.6 The lever on the righthand side of the Fine Adjust Orientating Head pulls forward to release the Cryocassette Clamp.

4.6.7 The knurled rotary control at the top of the Specimen Head sets the vertical orientation of the specimen. Clockwise rotation tilts the specimen upward; counter-clockwise rotation tilts the specimen downward.


4.6.8 The knurled rotary control at the lefthand side of the Specimen Head sets the horizontal orientation of the specimen. Clockwise rotation tilts the specimen to the left; counter-clockwise rotation tilts the specimen to the right.

4.6.9 Push the lock lever away from you (clockwise) to clamp the Specimen Head in the set position.


**Note**

- 1 *A SIMPLE ADJUST orientating head with clamp lever and joystick control is also available.*


4.6.4




4.6.5




4.6.6




4.6.7



4.6.8



4.6.9



# OPERATION


## 5.1 INTRODUCTION

5.1.1 When not in use, the instrument is normally kept switched on with the sliding window locked closed and the Refrigerated Chamber at approximately  $-15^{\circ}\text{C}$  to  $-20^{\circ}\text{C}$ .


5.1.2 The Control Panel is disabled, and the instrument is inoperable, when the keyswitch is vertical. The yellow **LOCKED** indicator is lit when the Control Panel is disabled and the key can be removed to keep the instrument locked.

5.1.3 Insert the key in the keyswitch and turn it clockwise to unlock the Control Panel. The **LOCKED** indicator is not lit when the Control Panel is enabled, and the key cannot be removed.


5.1.4 If [**WINDOW LOCK**] is lit, press the [**WINDOW LOCK**] push-button once to release the window and extinguish the indicator.



5.1.2



5.1.3



5.1.4

5.2.5 When the mountant starts to freeze, place the specimen directly from the cutting-up board on to the mountant and press down lightly to eliminate any entrapped air.


5.2.6 Add a final coat of mountant and allow it to freeze.

5.2.7 Lower the heat sink gently onto the specimen. Allow the specimen and mountant to freeze until both are opaque and firm. Take care not to touch the specimen.


#### Notes

- 1 Do not use too much mountant, or get it onto the sides or back of the Cryocassette.
- 2 Normal compressor operation resumes 10 minutes after cryobar boost is initially enabled. During this period, you may:
  - i manually switch off the cryobar boost facility by pressing **[CRYOBAR BOOST]**.
  - ii Enable Cryobar Boost for a further 10 minutes by pressing **[CRYOBAR BOOST]** while the Cryocar Boost LED is flashing.

5.2.6



5.2.7



## 5.4 TO CUT SECTIONS



**KEEP THE WINDOW CLOSED WHEN CUTTING OR INJURIOUS AEROSOLS MAY BE INHALED.**

5.4.1 At the Minor Control Panel use the **MICRON SELECTOR** and the [+] and [-] push-buttons to set the required section thickness.

5.4.2 Turn the Flywheel smoothly. The instrument performs a single cutting stroke with each full revolution of the Flywheel, and cuts a section of the specimen. Check to make sure that the section is flat before removing it onto a slide.


5.4.3 If necessary, adjust the Anti Roll Plate to ensure that each section is taken cleanly and is un-curved.

5.4.5 Continue to rotate the Flywheel if additional sections are required. The Specimen advances automatically by the set distance before the start of each cutting stroke.


## 5.5 TO REMOVE THE SPECIMEN

5.5.1 When the last section is cut, lock the Flywheel at the 3 o'clock position. Remove the cryocassette from the clamp. If the specimen is to be stored frozen, wrap it in metal foil before putting it into a freezer. **Do not store it in the Refrigerated Chamber of the Cryotome.**

5.4.1



5.4.3



# DECONTAMINATION CLEANING, AND MAINTENANCE


## 6.1 INTRODUCTION

6.1.1 It is important that decontamination and cleaning of the instrument becomes an automatic routine - especially if the source of material is unknown. It is recommended that a log is kept which lists:

- i Material Cut.
- ii Degree of Risk.
- iii Decontamination Performed.
- iv Name of User.
- v Department.

6.1.2 Cleaning is very much a part of good laboratory housekeeping practice. External cleaning requirements are straightforward. The Refrigerated Chamber should only be washed and cleaned after the Chamber is thoroughly decontaminated by fumigation.

6.1.3 Routine maintenance by the Operator does not require the removal of any panels or fixtures. Precision maintenance and adjustment should only be performed by Shandon trained Service personnel and it is recommended that a Maintenance Contract is obtained from our Service Department.



6.2.6 If DAY 9 is set in the fumigation program, the instrument goes into 'Standby Mode' (PSAVE) after fumigation is completed. Therefore, if a DAY 9 option is selected in the fumigation program, insert a beaker containing 10 ml of Ammonia SG880 in the Refrigeration Chamber to absorb residual formaldehyde from the atmosphere in the Chamber. Keep the chamber closed while allowing the beaker to stand for an hour.

### 6.3 CLEANING


6.3.1 Use a proprietary spray foam cleanser, or a cloth damped in soapy water, and a dry polishing cloth to periodically clean the outside of the cabinet. DO NOT USE EXCESSIVE WATER.

6.3.2 Periodically wash the Refrigerated Chamber with alcohol or warm soapy water to remove debris. Rinse well with a solution of 10% Formalin. Wash the Formalin solution down the outlet tubes in the bottom of the chamber, and beneath the fins of the heat exchanger.

6.3.3 Wipe and dry the walls of the Refrigerated Chamber carefully using a chamois leather damped with the cleaning solution followed by a dry polishing cloth.

6.3.4 The high polish finish inside the Refrigerated Chamber is designed to minimise the adhesion of particles to the sidewalls.

**WARNING** DO NOT USE SCOURING POWDER OR HARSH DETERGENTS OR YOU WILL DEGRADE THE POLISHED FINISH.



## 6.5 MAINTENANCE




**SWITCH OFF THE INSTRUMENT AT THE MAINS AND REMOVE THE PLUG FROM THE WALL SOCKET BEFORE PROCEEDING WITH ANY MAINTENANCE.**

### 6.5.1. TO REMOVE THE LAMP TUBE


- i Lift off the top cover,
- ii Loosen the two screws that hold the shade.
- iii Slide the shade along its two keyhole slots, then lift it off.
- iv Grip the fluorescent tube firmly, rotate it 90°, and pull it towards you to slide it from the end connector.

### 6.5.2 TO FIT A LAMP TUBE


- i Hold the replacement tube firmly at the ends and align the connector pins with the slots at each end of the fitting.
- ii Carefully push the tube into place.
- iii Rotate the tube a quarter turn.
- iv Fit the shade into its two keyhole slots and push it into position.
- v Tighten the two screws that hold the shade, and replace the cover.



6.5.1.i



6.5.1.ii



6.5.1.iv

# TROUBLE SHOOTING

7.1 Correct service and maintenance is essential for the long term serviceability of precision engineered products such as Cryotome®. We strongly recommend that a Shandon Service Contract is used to ensure future reliability, and consistency of performance.

7.2 Cryotome® is designed for maximum safety, ease of use, and reliability. Self test routines are built in so that fault codes (F:(Number)) show in the displays if malfunctions occur. No operator response is involved following a fault code other than to contact your Shandon Engineer for service. Other problems and their causes, that are not shown in the displays, are listed in Table 1. Alarm signals, their causes, and the remedial action to be taken are listed in Table 2.

7.3 If a fault indication continues or recurs, remove the samples from the refrigerated chamber and store them in a suitable refrigerated container.

**TABLE 1 - PROBLEMS WHEN SECTIONING**

SYMPTOM	CAUSE	REMEDY
Specimen cracks as it is frozen.	1 Too rapid freezing. 2 Specimen too thick.	1 Cut specimen into thinner block.
Specimen falls off Cryocassette.	1 Specimen too thick. 2 Cryocassette too cold before mountant added. 3 Specimen unsupported.	1 Cut specimen into thinner block. 2 Reduce time before adding mountant. (Should be at room temperature) 3 Give the specimen more support.
Specimen advances but does not cut.	1 Knife loose. 2 Cryocassette loose.  3 Specimen not firmly stuck to the Cryocassette. 4 Knife angle incorrect.	1 Make sure that the knife is clamped correctly 2 Make sure that the Cryocassette is fitted and clamped correctly. 3 Re-mount the specimen on the Cryocassette using adequate mountant. 4 Increase knife angle.
Sections roll-up.	1 Incorrect gap between Anti Roll Plate and knife. 2 Anti Roll Plate too low. 3 Incorrect Anti Roll Plate angle.	1 Adjust Anti Roll Plate gap.  2 Raise the height of the Anti Roll Plate. 3 Adjust the angle between the Anti Roll Plate and the knife.
Sections thaw when cut.	1 Cutting equipment not cold enough.	1 Allow more time for the knife and the Anti Roll Plate to cool.

**TABLE 2 - ERROR MESSAGES**

<b>SYMPTOM</b>	<b>CAUSE</b>	<b>REMEDY</b>
<b>LOBAT</b>	1 Internal back-up battery is losing charge	1 Switch on/leave the unit switched on. 2 Have the battery replaced.
<b>Err 1</b>	1 Invalid pushbutton entry. 2 Too many pushbuttons pressed simultaneously.	1 Use the correct pushbutton. 2 Re-enter pushbutton selections correctly.
<b>Err 2</b>	1 Invalid pushbutton entry during Defrost/ Fumigate. 2 Defrost is inhibited because the Cryobar and Chamber are too warm.	1 Use the correct pushbutton, or wait for the function to end. 2 Adjust working temperature to a lower chamber temperature.
<b>Err 3</b>	Not allocated on Electronic variants.	None
<b>Err 4</b>	1 Window lock error - cannot lock or unlock.	1 Make sure that the window is closed correctly. 2 Contact the Shandon Service Engineer.
<b>Err 5</b>	Not allocated on Electronic variants.	None.
<b>Err 6</b>	Not allocated on Electronic variants.	None.
<b>Err 7</b>	Not allocated on Electronic variants.	None
<b>Err8</b>	Not allocated on Electronic variants.	None

# SPECIFICATION

## Physical

Height(mm)	1070mm	(42 1/8 ins)
Width	660mm	(26 ins)
Depth	640mm	(25 1/4 ins)
Weight	150Kg	(330 lbs)
Mounting	4 off Movable castors - Front 2 retractable	

## Electrical

*NOTE Voltage is set at the factory.*

Power	50Hz = 1700 VA; 60Hz = 1400 VA
a.c. Voltage (V1)	220 - 240; 50 Hz
(V2)	110 - 120; 60 Hz
(V3)	220V 60 Hz

Switch convention

**I** = Power On  
**O** = Power Off

## Indicators

Display Panels  
Pushbuttons

3 off LCD - MICRONS, TIME and TEMPERATURE  
Illuminating - Yellow  
Steady State = Selected  
Flashing = Operating  
Specimen Travel Bar graph.

## Controls

Mechanical  
Electrical

Flywheel - Cutting Head drive  
Advance/Rewind Controls  
Micron Select  
Reset pushbutton  
Time Controls  
Chamber Controls

- Temperature control.
- Defrost.
- Fumigate.
- Lamp
- Cryobar
- Window Lock

Isolating keyswitch

## Capacities

Waste bottle	2 litres.
Fumigation Reservoir	4 ml

## Refrigeration

Conventional refrigeration system employing heat exchanger inside the refrigerated chamber.  
Minimum Chamber Temperature - 30°C at 20°C ambient or lower.

## **Basic Equipment**

DESCRIPTION		PART NUMBER
Cryotome E	220/240 V 50 Hz	0620E
	110/120 V 60 Hz	0620E/110
	220 V 60 Hz	0620E/220
Plastic Mat		0620-192
2 litre bottle		P13198

## **Accessories Pack supplied with Basic Equipment**

DESCRIPTION	PACK OF:	PART NUMBER
Cryocassettes	5	0620-001
Shelf		0620-002
Debris Tray		0620-004
Flywheel Assembly		0620-118
Flywheel Bolt		0300-095
Flywheel Bolt Key		0300-124
Mains Lead	U.K	P13291
	EU	P13290
	U.S	P13292
Brush		P12940
Front Panel Keys	2	P13571
Operator Guide		77210163

## **Accessories Available on Separate Order**

Solid Knife holder with Anti Roll Plate	0620-023
High Profile Disposable Knife holder with Anti Roll Plate	0620-021H
Low Profile Disposable Knife holder with Anti Roll Plate	0620-021L

## **Cryocassettes**

Cryocassettes	5	0620-001
Biopsy Cryocassettes	5	0620-007
Heavy Duty Cryocassettes	5	0603
Specimen Tube	5	0605
Large Specimen Holder Tube (50mm)		0605-050
Simple Orientating Head		0620-008
Fine Adjust Orientating Head		0620-006
Specimen Tube Adaptor		0620-820
Cryocassettes (28x40mm)	5	0620-026
(45mm dia)	5	0620-036
(35x35mm)	5	0620-039
(35mm dia)	5	0620-041
(Assorted)	4	77210080

Continued

# WARRANTY STATEMENT

9.1.1 Provided the terms of payment are duly complied with the Company undertakes to remedy any original defects arising from faulty materials or workmanship in any goods manufactured by the Company which under proper and normal conditions of use may develop within a period of twelve months from the date of delivery provided the same are returned to the Company as provided by this paragraph.

9.1.2 In the case of components which by their nature of application have an unpredictable life, this guarantee shall only be to the extent of the guarantee given by the manufacturers of these articles.

9.1.3 The Company will accept no liability where in the opinion of the Company the defect has been caused by damage due to the Buyer's failure to follow operating instructions, incorrect installation, wear and tear, or damage due to the use of spare parts other than those spare parts of the Company or which are recommended by the Company or where in the opinion of the Company the defect has been caused by alterations or repairs being undertaken by any person other than the Company.

9.1.4 Any damage claim must be in writing and give the serial number and description of goods, order number and date of delivery and will not apply where any names or serial numbers or other information which may be attached to or inscribed upon the goods have been removed, covered up or defaced in any way.

9.1.5 Any goods or parts thereof which may require repair or replacement shall be repaired or replaced (at the election of the Company) at the works of the Company only to which the Buyer shall deliver the same carriage paid at their risk and at the Buyer's expense. Any such goods or parts will be delivered by the Company to the Buyer free within the United Kingdom but if required to be delivered elsewhere the freight insurance and other charges from works to destination shall be borne by the Buyer. All faulty parts removed from the equipment will become the Company's property. Any other repairs or work by the Company will be carried out under the terms and conditions for specialist engineers currently in force.

9.1.6 If any goods or parts thereof are returned unnecessarily all costs involved, including a charge for inspection, handling and the return carriage must be paid by the sender. In no circumstances shall any of the goods be returned to the Company without its prior written consent.

# INDEX

	TOPIC	SECTION	PAGE
<b>A</b>	Advance .....	1.1.14 .....	1.4
	[ADVANCE] .....	1.1.20 .....	1.6
	ALd .....	3.11.1.ii .....	3.11
	Anti-Roll Plate .....	1.1.9 .....	1.3
	Anti Roll Plate, To Set Up the Knife Holder and ~ ..	4.5 .....	4.6
<b>B</b>			
<b>C</b>			
	Cancel Defrost/Fumigate .....	3.8.8 .....	3.8
	Carriage .....	1.1.11 .....	1.3
	Chamber Controls .....	3.8 .....	3.6
	Cleaning .....	6.3 .....	6.3
	<b>CLEANING AND MAINTENANCE, DECONTAMINATION ~ .....</b>		<b>6.1</b>
	Cleaning .....	6.3 .....	6.3
	Defrosting .....	6.4 .....	6.4
	Fumigation .....	6.2 .....	6.2
	Maintenance .....	6.5 .....	6.5
	To Fit a Lamp Tube .....	6.5.2 .....	6.5
	To Remove a Lamp Tube .....	6.5.1 .....	6.5
	Clock, To Set The ~ .....	3.5 .....	3.4
	Control Panel, Main ~ .....	1.1.19 .....	1.6
		3.1.1 .....	3.1
	<b>CONTROLS.....</b>		<b>3.1</b>
	Chamber Controls .....	3.8 .....	3.6
	Cancel/Defrost/Fumigate .....	3.8.4 .....	3.6
	Cryobar Boost .....	3.8.10 .....	3.8
	Immed Defrost .....	3.8.3 .....	3.6
	Immed Fumigate .....	3.8.6 .....	3.7
	Lamp .....	3.8.9 .....	3.8
	Window Lock .....	3.8.12 .....	3.8
	Description.....	3.1 .....	3.1
	Main Control Panel .....	3.1.1 .....	3.1
	Minor Control Panel .....	3.1.3 .....	3.1
	Micron Selector .....	3.1.3 .....	3.1
	Rewind/Advance .....	3.1.4 .....	3.1
	Keyswitch Control .....	3.10 .....	3.10
	Message Codes .....	3.11 .....	3.10
	Lobat .....	3.11.1.i .....	3.10
	Ald .....	3.11.1.ii .....	3.11
	dEF .....	3.11.1.iii .....	3.11
	--- .....	3.11.1.iv.....	3.11

Continued

Continuation

	<b>TOPIC</b>	<b>SECTION</b>	<b>PAGE</b>
<b>F</b>	Fault Codes (F:)	3.11.1.v	3.11
	Flywheel	1.1.12	1.3
	Fumigation	6.2	6.2
	Fumigate, To Set ~ Times	3.7	3.5
<b>G</b>			
<b>H</b>	Housekeeping	5.6	5.6
	Daily	5.6.1	5.6
	Daily after Fumigation	5.6.2	5.6
<b>I</b>	Immed Defrost	3.8.3	3.6
	[IMMED FUMIGATE]	1.1.6	1.5
		3.8.6	3.7
	<b>INDEX</b>		<b>10.1</b>
	<b>INSTALLATION INSTRUCTIONS</b>		<b>2.1</b>
	Procedure	2.1	2.1
<b>J</b>			
<b>K</b>	Keyswitch Control	3.10	3.10
	Locked Indicator	3.10.2	3.10
	Knife Holder	1.1.8	1.2
	Knife Holder, To Fit a ~	4.2	4.2
	Knife Holder, To Set the ~ and Anti-Roll Plate	4.5	4.6
	Knife Holder, To Set the ~ for Cutting Specimen	5.3	5.4
<b>L</b>	[LAMP]	1.1.21	1.6
		3.8.9	3.8
	Lamp, To Fit a ~ Tube	6.5.2	6.5
	To Remove the ~ Tube	6.5.1	6.5
	Lobat	3.11.1.i	3.10
	'Locked' Indicator	3.9.1	3.9
	5.1.3	5.1	
<b>M</b>	Main Control Panel	3.1.1	3.1
	Maintenance	6.5	6.5
	To Remove the Lamp Tube	6.5.1	6.5
	To Fit the Lamp Tube	6.5.2	6.5

Continued

Continuation

	TOPIC	SECTION	PAGE
<b>Q</b>			
<b>R</b>	Remove the Specimen, To .....	5.5 .....	5.6
	Reset .....	1.1.15 .....	1.4
	rESEt .....	3.11.2.ii .....	3.11
	Retracted .....	1.1.13 .....	1.4
	Rewind/Advance .....	1.1.20 .....	1.6
<b>S</b>	Safety.....	3 .....	3
	Sections, To Cut ~ .....	5.4 .....	5.5
	<b>SETTING UP</b> .....		<b>4.1</b>
	Introduction .....	4.1.1 .....	4.1
	Specimen Head Controls .....	4.6 .....	4.8
	To Fit a Disposable Blade .....	4.4 .....	4.4
	To Fit a Knife Holder .....	4.2 .....	4.2
	To Fit a Standard Solid Knife .....	4.3 .....	4.3
	To Set Up the Knife Holder and Anti Roll Plate .....	4.5 .....	4.6
	Set Clock.....	3.5 .....	3.4
	Set Chamber Temperature .....	3.9 .....	3.9
	Solid Knife, To Fit a ~ .....	4.3 .....	4.3
	Spares List .....		8.3
	<b>SPECIFICATION</b> .....		<b>8.1</b>
	Accessories .....		8.3
	Basic Equipment .....		8.3
	Specimen Head .....	1.1.4 .....	1.1
		1.1.7 .....	1.2
	Specimen Head Controls .....	4.6 .....	4.8
	Orientation, Cryocassette ~ Clamp .....	4.6.4 .....	4.9
	Cryocassette Clamp .....	1.1.7 .....	1.2
	Specimen, Mounting ~ .....	5.2 .....	5.2
	Specimen, To Remove ~ .....	5.5 .....	5.6
	Specimen Travel Available .....	3.3 .....	3.3
	STANDBY (PSAVE), To Set .....	3.6.3 .....	3.5
<b>T</b>	Table 1 - Problems When Sectioning .....		7.1
	Table 2 - Error Messages .....		7.3
	TEMPERATURE ° C Display .....	3.9.1 .....	3.9
		3.11.1 .....	3.10
	Temperature Controls .....	3.9 .....	3.9

Continued

## Declaration of Conformity

**Manufacturer's Name:** Life Sciences International (Europe) Limited,  
**Manufacturer's Address:** 93 - 96 Chadwick Road, Astmoor, Runcorn,  
Cheshire WA7 1PR  
ENGLAND  
**Product Description:** Electronic Cryostat  
**Product Designation:** **Cryotome® E**  
**Year of Marking (CE):** 1996

This product conforms with the essential protection requirements of the EMC Directive,  
89/336/EEC (*as amended by 92/31/EEC & 93/68/EEC*),  
and the Low Voltage Directive 73/23/EEC (*as amended by 93/68/EEC*).

This product complies with the following International Standards:

**EMC:** EN55022 class A  
EN50082 - 1

**Safety** IEC1010 - 1, 1990 +Amd.1,1992 + Amd.2,1995

**Issued by:** Mr R.Russell-Smith  
Quality Assurance Manager  
Life Sciences International (Europe) Limited



(Signature)